

CLEANING AND CARE OF CENTRIFUGES



Regular cleaning and maintenance ensures the long lifespan of your Hettich centrifuge. Reduce unnecessary repair costs by following these recommendations.

We suggest setting a regular time daily or weekly – only 15 minutes per week ensure that the failure probability of your centrifuge is significantly reduced. This poster does not replace the operating manual.

PREVENTIVE ACTIONS



Step 1
Switch off the device and disconnect the centrifuge from the power supply.

Step 2
For cleaning, use appropriate PPE (personal protective equipment, e.g. rubber gloves).

Step 3
All you need for cleaning is a mild detergent with a pH value of 5-8 or a 70% alcohol solution and slightly damp paper towels.

Step 4
Clean the surface of the device (housing).

