HettInfo II

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Andreas Hettich GmbH & Co. KG
Föhrenstraße 12, D-78532 Tuttlingen / Germany
Phone +49 (0)7461 / 705-0
Fax +49 (0)7461 / 705-1125
info@hettichlab.com, service@hettichlab.com
www.hettichlab.com

Single Registration Number:

DE-MF-000010680

© 2020 by Andreas Hettich GmbH & Co. KG

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
OperatorID.txt

11 Free and Open Source Software

This product contains free and open source software. See the enclosed DVD "HettInfo II Licenses and Sourcecode" for the applicable license texts, the source code, copyright notices for the individual components, any disclaimers, permission notes and other information about the individual components..

The DVD contains the following files/folders:

File	Content
File u-boot-*.*	Folders with the license text of the uboot bootloader components
gelin2-21.10-ix6-x86_64.tar.*	Complete Linux packet (GELinBSP) incl. license info and source code
HettInfo-GELin_KernelV5.4.196_Lizenzinformationen.txt	Links to the license info of the Linux Kernel V5.4.196
hettinfo-gelin_21.10-bom.csv	List of all GELin components with detailed information to the licenses
hettinfo-gelin-license-information.*	License text files (txt and HTML)
linux-5.4.196-ge-22.04+-exceed_ts-hmi-hettich.zip	Linux Kernel V5.4.196 (adapted for the use with the HMI electronic board)

Please note: Use of free and open source software components listed on the DVD files mentioned in the table above are governed by the terms of the licenses of the respective copyright holder as described in the folders on the DVD. Some of the licenses contain liability and warranty disclaimers and are provided "as is" without warranty of any kind. For more details and additional conditions for using the software, please read the attached licenses.

License Terms and Conditions regarding LGPL components:

For those parts of the HETTICH software that are licenses under

- the "GNU Library General Public License" or the "GNU Lesser General Public License" (LGPL) version 2.1 and earlier ("work") the customer is expressly allowed to make modifications of the work for the customer's own use and to reverse engineer for debugging such modifications.
- the "GNU Lesser General Public License" (LGPL) version 3 and later ("work") it is expressly allowed to make modifications of the portions of the library contained in the work and reverse engineer for debugging such modifications.

Further rights of the customer under applicable copyright law are not restricted. Please note, that those parts of the software which are under free and open source licenses are distributed free of charge.

Installation Information for modified components:

Modified software components, esp. those under the GNU Lesser General Public License version (LGPL), can be installed by copying the new software component to the storage device contained within the product.

	CAUTION: CHANGING THE SOFTWARE OF THE PRODUCT MAY RESULT IN INJURY TO LIFE, LIMB, OR HEALTH! ALSO IT CAN IMPACT THE CERTIFICATION OF THE PRODUCT AND RESULT IN A NONLEGAL USE OF THE DEVICE. CHANGES CAN ALSO AFFECT THE WARRANTY.
---	---

Special License Information

The enclosed DVD "HettInfo II Licenses and Sourcecode" is part of the present documentation. We would like to point out the following separately:

OpenSSL

This product includes software developed by the OpenSSL Project for use in the OpenSSL Toolkit (<http://www.openssl.org>).

Eric Young

This product includes cryptographic software written by Eric Young (ey@cryptsoft.com).

Fontconfig

This product includes software under the fontconfig license.

- Copyright © 2000,2001,2002,2003,2004,2006,2007 Keith Packard
- Copyright © 2005 Patrick Lam
- Copyright © 2009 Roozbeh Pournader

-
- Copyright © 2008,2009 Red Hat, Inc.
 - Copyright © 2008 Danilo Šegan
 - Copyright © 2012 Google, Inc.

Permission to use, copy, modify, distribute, and sell this software and its documentation for any purpose is hereby granted without fee, provided that the above copyright notice appear in all copies and that both that copyright notice and this permission notice appear in supporting documentation, and that the name of the author(s) not be used in advertising or publicity pertaining to distribution of the software without specific, written prior permission. The authors make no representations about the suitability of this software for any purpose. It is provided "as is" without express or implied warranty.

THE AUTHOR(S) DISCLAIMS ALL WARRANTIES WITH REGARD TO THIS SOFTWARE, INCLUDING ALL IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS, IN NO EVENT SHALL THE AUTHOR(S) BE LIABLE FOR ANY SPECIAL, INDIRECT OR CONSEQUENTIAL DAMAGES OR ANY DAMAGES WHATSOEVER RESULTING FROM LOSS OF USE, DATA OR PROFITS, WHETHER IN AN ACTION OF CONTRACT, NEGLIGENCE OR OTHER TORTIOUS ACTION, ARISING OUT OF OR IN CONNECTION WITH THE USE OR PERFORMANCE OF THIS SOFTWARE.

The FreeType Project License

This product includes software based in part of the work of the FreeType Team (<https://freetype.org/>).
For further information (including license texts and copyright notices) see the enclosed DVD.

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All rights reserved. No part of this publication may be reproduced without the prior written permission of the copyright owner.