

ROTIXA 500 RS



Inhalt des Dokuments / content of the document

Operating instructions (EN)

Rotoren und Zubehör / Rotors and accessories

Operating instructions

ROTIXA 500 RS



Translation of the original operating instructions

©2023 – All rights reserved

Andreas Hettich GmbH & Co. KG

Föhrenstrasse 12

D-78532 Tuttlingen, Germany

Telephone: +49 (0)7461 705-0

Fax: +49 (0)7461 705-1125

Email: info@hettichlab.com, service@hettichlab.com

Internet: www.hettichlab.com

Table of contents

1	About this document	6
1.1	Use of this document	6
1.2	Gender reference	6
1.3	Symbols and labels in this document	6
2	Safety	6
2.1	Intended use	6
2.2	Personnel requirements	7
2.3	Operator's responsibility	8
2.4	Safety instructions	8
3	Device overview	10
3.1	Technical data	10
3.2	European registration	13
3.3	Important labels on the packaging	14
3.4	Important labels on the device	14
3.5	Operating and indicator elements	16
3.5.1	Control	16
3.5.2	Indicator elements	16
3.5.3	Controls	17
3.6	Original spare parts	19
3.7	Scope of supply	19
3.8	Returns	19
4	Transport and storage	20
4.1	Transport and storage conditions	20
5	Commissioning	21
5.1	Unpacking the centrifuge	21
5.2	Setting up and connecting the centrifuge	22
5.3	Switching the centrifuge on and off	24
6	Operation	24
6.1	Opening and closing the lid	24
6.2	Removing and installing the rotor	25
6.3	Inserting and removing buckets	26
6.4	Inserting and removing adapters	26
6.5	Loading	27
6.6	Opening and closing the biosafety system	29
6.6.1	Explanation	29
6.6.2	Lid with screw cap and hole	30
6.6.3	Lid with screw cap	30
6.7	Packing instructions, HettLiner	30
6.8	Centrifugation	32
6.8.1	Centrifugation in continuous operation	32
6.8.2	Centrifugation with time preselection	33
6.8.3	Changing settings during centrifugation	33
6.9	Quick stop function	33

7	Software operation.	33
7.1	Key switches.	33
7.2	Centrifugation parameters.	34
7.2.1	Ramp-up and ramp-down parameters.	34
7.2.2	Runtime.	34
7.2.3	Speed, RPM.	35
7.2.4	Integral RCF.	35
7.2.5	temperature	35
7.2.6	Relative centrifugal force, RCF.	35
7.2.7	Setting the relative centrifugal force (RCF/RZB).	36
7.2.8	Centrifugation of substances or mixtures of substances with a density higher than 1.2 kg/dm ³ .	36
7.2.9	Centrifuging radius.	36
7.3	Programming.	36
7.3.1	Opening or loading programs.	36
7.3.2	Entering or changing programs.	37
7.3.3	Automatic buffer.	37
7.4	Rotor detection.	37
7.5	Cooling (for centrifuges with cooling).	37
7.5.1	Instructions, cooling.	37
7.5.2	Standby cooling.	37
7.5.3	Precooling the rotor.	37
7.6	Heating (for centrifuges with heating).	38
7.7	Machine Menu.	38
7.7.1	Querying system information.	38
7.7.2	Querying operating hours.	39
7.7.3	Audible signal.	39
7.7.3.1	General.	39
7.7.3.2	Setting an audible signal.	39
7.7.4	Centrifugation data displayed after switching on.	39
7.7.5	Setting the date and time.	40
7.8	Program links.	40
7.8.1	Linking programs or changing a program link.	40
7.8.2	Centrifugation run with program link.	41
7.8.3	Deleting program links.	42
8	Cleaning and care.	42
8.1	Overview table.	42
8.2	Cleaning and disinfection instructions.	43
8.3	Cleaning.	43
8.4	Disinfection.	44
8.5	Maintenance.	45

9	Troubleshooting	46
9.1	Fault description	46
9.2	Perform a MAINS RESET	47
9.3	Emergency release	48
10	Disposal	48
10.1	General instructions	48
11	Index	50

1 About this document

1.1 Use of this document

- Read this document carefully and in full before commissioning the device for the first time.
Observe other enclosed instruction sheets where necessary.
- This document is part of the device and must be kept within easy reach.
- This document must be included if the device is passed on to a third party.
- The current version of the document in the available languages can be found on the manufacturer's website: ➔ <https://www.hettichlab.com/de/download-center/>

1.2 Gender reference

The employed masculine or feminine language form is to facilitate reading. In the spirit of equal treatment, corresponding terms apply in principle to all genders and do not imply any valuation.

1.3 Symbols and labels in this document

General symbols

The following markers are used in this document to highlight instructions, results, listings, references and other elements:

Marker	Explanation
1. → 2. → 3. → ... →	Step-by-step instructions
➡	Results of action steps
➡	References to sections of the document and other applicable documents
■ ... ■ ...	Listings without a fixed order
/Buttons/	Controls (for example: buttons, switches)
'Indicator'	Indicator elements (for example: signal lights, screen elements)

2 Safety

2.1 Intended use

Intended use

This device is a laboratory centrifuge suitable for medical applications.

Their exclusive therapeutic purpose is to centrifuge blood in blood bag systems. The separated blood components are transferred by another device (separator) into corresponding satellite bags. The individual components obtained in this way are then used for transfusion or autotransfusion.

The centrifuge is only to be operated by qualified personnel working for blood donation services or hospitals.

The centrifuge is only intended for the uses referred to above.

Any other use or use beyond this is considered improper. Andreas Hettich GmbH & Co. KG shall not be liable for any damage arising from this.

Intended use also includes the observation of all instructions in the user manual and compliance with the required inspection and maintenance intervals.

Non-intended use

- The centrifuge is not suitable for use in explosive or radioactive, or biologically or chemically-contaminated atmospheres.
- The user must take appropriate actions when centrifuging hazardous substances or mixtures of substances that are toxic, radioactive or contaminated with pathogenic microorganisms.
The manufacturer generally recommends using only centrifuge tubes with special screw caps for hazardous substances.
Use sealable centrifuge tubes with a biosafety system for materials of risk groups 3 and 4.
- The manufacturer does not recommend centrifugation with flammable or explosive materials.
- The manufacturer does not recommend centrifugation with materials that react chemically with one another with high energy.

Foreseeable misuse

The manufacturer recommends using only accessories that it has approved for the intended purpose.

Only operate the centrifuge under supervision.

2.2 Personnel requirements

Required qualifications

The user has read the user manual in full and familiarised themselves with the device.



NOTICE

Damage to the device by unauthorised personnel

- Tampering with and modifications to devices by unauthorised persons are at the operating organisation's own risk and will result in the loss of all warranty and liability claims.

Trained user

The user is trained in laboratory practice and able to carry out the work assigned to them, and to recognise and prevent potential hazards independently.

Personal protective equipment

Lack of personal protective equipment or unsuitable personal protective equipment increases the risk of impaired health and injury.

- Only use personal protective equipment that is in proper condition.
- Only use personal protective equipment that is adapted to the person (correct size, for example).
- Observe instructions on other protective equipment for specific activities.

2.3 Operator's responsibility



Follow the instructions in this document for proper and safe use of the device.

Keep the user manual for future reference.

Provide information

- Following the instructions in this document will help:
 - To avoid dangerous situations.
 - To minimise repair costs and downtime.
 - To increase the reliability and service life of the device.
- The operator is responsible for compliance with company regulations, standards and national laws.
- Note and keep the revision of the document separate from the document. If lost, the document can be replaced in the correct revision.
- Keep the user manual available at the place where the device is used.
- Pass the user manual on to the buyer when the device is sold.

Training of personnel

Lack of knowledge when working with the device may result in serious injury or death.

- Instruct personnel on their tasks and the associated risks in accordance with the instruction.

2.4 Safety instructions



Reporting serious incidents and notifiable incidents

In the event of serious incidents or notifiable incidents involving the device or its accessories, these must be reported to the manufacturer and, where applicable, to the competent authority where the user and/or the patient is registered.



DANGER

Risk of contamination for the user due to inadequate cleaning or failure to observe the cleaning instructions.

- Observe cleaning instructions.
- Wear personal protective equipment when cleaning the device.
- Observe laboratory regulations (e.g. TRBAs, the German Protection against Infection Act, hygiene plan) for handling biological agents.



DANGER

Fire and explosion hazard due to hazardous substances in samples.

- Observe relevant regulations and directives for handling chemicals and hazardous substances.
- Do not use aggressive chemicals (for example: dangerous, corrosive extraction agents such as chloroform, strong acids).

**WARNING**

Dangers due to insufficient maintenance or maintenance not carried out on time.

- Follow maintenance intervals.
- Check the device for visible damage or defects.
If any visible damage or defects are present, take the device out of service and inform a service technician.

**WARNING**

Risk of electric shock due to ingress of water or other liquids.

- Protect the device against external liquids.
- Do not pour any liquids into the interior of the device.
- Transport using original transport packaging.

**WARNING**

Contamination with hazardous substances and substance mixtures!

Observe the following actions for substances and substance mixtures that are toxic, radioactive and/or contaminated with pathogenic microorganisms:

- As a rule, use only centrifuge tubes with special screw caps for hazardous substances.
- Use sealable centrifuge tubes with a biosafety system for materials of risk groups 3 and 4.
- If no biosafety system is used, the device is not microbiologically tight in the sense of standard EN / IEC 61010-2-020.
- Contact the manufacturer if necessary.

**WARNING**

Risk of injury and damage to the device due to a loose rotor.

- The driver of the rotor shaft must be correctly seated in the groove of the rotor when mounting the rotor.
- Hand-tighten the nut securing the rotor.
- Check that the rotor is firmly seated.
- Follow maintenance intervals.

**CAUTION**

Risk of injury due to rotating rotor

Long hair and items of clothing can get caught on the rotor if the rotor is moved manually.

- Tie long hair back.
- Do not allow garments to hang in the centrifuging chamber.

**NOTICE**

Damage to the device electronics due to incorrect voltage or frequency at the device circuit breaker.

- Operate the device with the correct mains voltage and mains frequency.

The value can be found in the technical data and on the rating plate.

**NOTICE**

Damage to the device and samples due to premature program termination.

Premature program termination is caused by power failure, switching off during the program or pulling out the mains plug.

- Do not switch off the device while the program is running.
- Do not trigger the emergency release on the device while the program is running.
- Do not pull out the mains plug while the program is running.

3 Device overview

3.1 Technical data

Manufacturer	Andreas Hettich GmbH & Co. KG, D-78532 Tuttlingen			
Model	ROTIXA 500 RS			
Type	4950, 4950-50		4950-70, 4950-80	
Mains voltage ($\pm 10\%$)	230-240 V 1~	220 V 1~	230-240 V 1~	220 V 1~
Mains frequency	50 Hz	60 Hz	50 Hz	60 Hz
power consumption	3800 VA			
Power consumption	16 A			
Refrigerant	R452A			
max. capacity	4 x 1000 ml			
max. permissible density	1.2 kg/dm ³			
max. speed (RPM)	11500			
max. acceleration (RCF)	18038			
max. kinetic energy	59620 Nm			

Obligation to perform checks (DGUV Rules 100-500) (valid only in Germany)	yes
---	-----

Ambient conditions (EN / IEC 61010-1):

Installation site	indoors only
Altitude	up to 2000 m above sea level
Ambient temperature	5 °C to 35 °C
Humidity	maximum relative humidity 80% for temperatures up to 31 °C, decreasing linearly to 50% relative humidity at 40 °C.
Overvoltage category (IEC 60364-4-443)	II
Pollution level	2
Device protection class	I not suitable for use in potentially explosive atmospheres.

EMC:

Emitted EM interference, EM interference immunity	EN / IEC 61326-1 Class B
Noise level (rotor-dependent)	≤65 dB(A) ≤63 dB(A)

Dimensions:

Width	650 mm
Depth	814 mm
Altitude	973 mm
Weight	approx. 219 kg approx. 233 kg

Manufacturer	Andreas Hettich GmbH & Co. KG, D-78532 Tuttlingen	
Model	ROTIXA 500 RS	
Type	4950-08, 4950-58 4950-78, 4950-88	
Mains voltage ($\pm 10\%$)	208 V 1~	
Mains frequency	60 Hz	
power consumption	3800 VA	
Power consumption	18 A	

Refrigerant	R452A
max. capacity	4 x 1000 ml
max. permissible density	1.2 kg/dm ³
max. speed (RPM)	11500
max. acceleration (RCF)	18038
max. kinetic energy	59620 Nm
Obligation to perform checks (DGUV Rules 100-500) (valid only in Germany)	yes

Ambient conditions (EN / IEC 61010-1):

Installation site	indoors only
Altitude	up to 2000 m above sea level
Ambient temperature	5 °C to 35 °C
Humidity	maximum relative humidity 80% for temperatures up to 31 °C, decreasing linearly to 50% relative humidity at 40 °C.
Overvoltage category (IEC 60364-4-443)	II
Pollution level	2
Device protection class	I not suitable for use in potentially explosive atmospheres.

EMC:

Emitted EM interference, EM interference immunity	FCC Class B	
Noise level (rotor-dependent)	≤65 dB(A)	≤63 dB(A)

Dimensions:

Width	650 mm	
Depth	814 mm	
Altitude	973 mm	
Weight	approx. 225 kg	approx. 239 kg

Rating plate




Fig. 1: Rating plate

- 1 Item number
- 2 Serial number
- 3 Revision
- 4 Equipment number
- 5 Data matrix code
- 6 any labelling indicating whether medical device or in vitro diagnostic medical device
- 7 Global Trade Item Number (GTIN)
- 8 Date of manufacture
- 9 Serial number
- 10 any EAC mark, CE mark
- 11 Country of manufacture
- 12 Date of manufacture
- 13 Mains frequency
- 14 Maximum kinetic energy
- 15 Maximum permissible density
- 16 Manufacturer's address
- 17 any Coolant circuit pressure
- 18 any Coolant capacity
- 19 any Coolant type
- 20 Revs per minute
- 21 Performance values
- 22 Mains voltage
- 23 any Device designation
- 24 Manufacturer's logo

3.2 European registration

Device conformity



Device conformity according to EU directives.

Notified body:

mdc medical device certification GmbH – Notified Body CE 0483

Tel: +49 (0)711 253597 0

Fax: +49 (0)711 258597 10

E-mail: mdc@mdc-ce.de

Website: www.mdc-ce.de

Address: Kriegerstrasse 6, D-70191 Stuttgart, Germany

Single Registration Number

SRN: DE-MF-000010680


Basic-UDI-DI

Basic-UDI-DI

Device assignment

040506740100029L

ROTIXA 500 RS (medical device)

3.3 Important labels on the packaging**TOP**


This is the correct upright position of the shipping container for transport and/or storage.

**FRAGILE GOODS**

The contents of the shipping container are fragile, so it must be handled with care.

**PROTECT FROM MOISTURE**

The shipping container must be kept away from rain and kept in dry conditions.

**TEMPERATURE LIMITATION**

The shipping container must be stored, transported and handled within the indicated temperature range (-20 °C to +60 °C).

**HUMIDITY LIMITATION**

The shipping container must be stored, transported and handled within the indicated air humidity range (10% to 80%).

**STACK LIMITATION BASED ON QUANTITY**

Maximum number of identical packages that may be stacked on the lowest package, "n" standing for the number of packages allowed. The lowest package is not included in "n".

3.4 Important labels on the device

The labels on the device must not be removed or covered, or have anything pasted over them.



Attention, general danger area.

Ensure you read the instructions for commissioning and operation and observe the safety instructions before using the device.



Biohazard warning.



Warning: hot surface.

Failure to observe this warning may result in damage to property and/or personal injury.



Direction of rotation of the rotor.

The orientation of the arrow indicates the rotor's direction of rotation.



Symbol for the separate collection of electrical and electronic equipment, in accordance with Directive 2012/19/EU (WEEE).

Use in European Union countries, Norway and Switzerland.



Key switch positions.



The centrifuge is equipped with an optical interface.

The optical interface is marked with a symbol.

The centrifuge can be controlled and data retrieved via the interface. The *[PROG]* button lights up during data communication.



Equipotential: Connector (PE plug) for equipotential bonding (only for centrifuges with a PE plug).

3.5 Operating and indicator elements

3.5.1 Control




Fig. 2: Control

3.5.2 Indicator elements

- The button lights up when the lid is closed.



Fig. 3: [Lid] button

- The indicator appears when the lid is closed.



Fig. 4: 'Lid closed' indicator

- The indicator appears when the lid is open.



Fig. 5: 'Lid open' indicator

LOCK 1, LOCK 2

Fig. 6: [Key switch position]
indicator

- The indicator appears when the key switch is in this switch position.

LOCK 4, LOCK 5

Fig. 7: [Key switch position]
indicator

- The indicator appears when the program lock is enabled during serial communication (for centrifuges with serial communication only).

PC,

Fig. 8: [Serial communication]
indicator

- The indicator appears if the centrifuge has a serial interface and the centrifuge is connected respectively not connected.



Fig. 9: 'Rotation' indicator

- The indicator appears when the rotor is turning.

STOP

Fig. 10: [STOP] indicator

- The indicator appears during the centrifugation run as long as the rotor is turning.
- The indicator flashes after an emergency stop.
- The indicator flashes after an emergency stop.

3.5.3 Controls



Fig. 11: [Rotary knob]

- Setting the individual parameters.
Turning anticlockwise decreases the value.
Turning clockwise increases the value.



Fig. 12: [Mains switch]

- Switch the device on and off.




Fig. 13: [Key switch]

- The key switch switches various functions on and off, depending on the position.




Fig. 14: [Temperature and centrifuging radius] button

- Temperature setpoint, parameter T/°C
Adjustable from -20°C to +40°C, in 1°C increments (adjustable from -20°C to +60°C with heating/cooling option).
The lowest achievable temperature is rotor dependent.
- Centrifuging radius
Parameter r/mm. Input in mm.



Fig. 15: [Ramp-up parameters] button

- Ramp-up levels, parameters
Level 9 = shortest ramp-up time, Level 1 = longest ramp-up time.
- Ramp-up time parameters
The adjustable time range is dependent on the set speed.



Fig. 16: [Ramp-down parameters] button

- Brake levels, parameters
R = Linear braking curve,
B = similar to an exponential braking curve.
Level R9, B9 = short ramp-down time, ...
Level R1, B1 = long ramp-down time,
Level R0 = unbraked ramp-down.
- Ramp-down time, parameters
The adjustable time range is dependent on the set speed.
- Brake cut-off speed, parameter n^(•) /RPM
Unbraked ramp-down takes place after reaching this speed.




- Open the lid.

Fig. 17: [Lid] button



- Querying the integral RCF, parameter fRCF

Fig. 18: [fRCF] button



- Speed, parameter RPM.
Adjustable from 50 RPM to the maximum rotor speed (n-max-Rotor)
- Querying the maximum rotor speed, parameter n-max-Rotor

Fig. 19: [n] button



- Select program location, parameter PROG No.

Fig. 20: [PROG] button



- Relative centrifugal force, parameter RCF/RZB.
A numerical value can be set that gives a speed between 50 RPM and the maximum rotor speed (n-max-Rotor). Adjustable in 1 second increments.
- Querying the maximum RCF of the rotor, parameter RCF-max-Rotor.

Fig. 21: [RCF] button



- Retrieving programs.

Fig. 22: [RCL] button



- Start centrifugation run.
- Acceptance of changes during the centrifugation run.

Fig. 23: [START] button



- Saving programs. 89 programs can be saved (program locations 1 to 89).
The program locations "----" and 90 to 99 serve as an automatic buffer.
No programs can be stored in these program locations.

Fig. 24: [STO] button



Fig. 25: [STOP] button



Fig. 26: [t] button

3.6 Original spare parts

Only use original spare parts from the manufacturer and approved accessories.

3.7 Scope of supply

The following accessories are supplied with the centrifuge:

- 1 grease for the trunnions
- 1 single-end open-ended spanner (10 mm AF)
- 1 open-ended spanner (17 mm and 19 mm AF)
- 1 Hex key (5 mm x 170)
- 1 square spanner
- 10 covering caps Ø12

- 3 wood screws
- 3 washers
- 2 metal rails
- 4 wide head nails
- 1 power cable
- 1 user manual
- 3 program data sheets for S control unit

Additionally for types 4950-70, 4950-78, 4950-80 and 4950-88:

- 1 Notes on setup and installation

Additionally for delivery in Germany:

- 1 inspection book

Rotors and the corresponding accessories are supplied depending on the order.

3.8 Returns

An original Return Material Authorisation (RMA) form from the manufacturer must always be requested for a return. Secure and reliable acceptance and booking in of the goods with the manufacturer is not possible without an

original RMA form from the manufacturer. The Return Material Authorisation (RMA) form contains a Declaration of No Objection (UBE), which must be completed in full and enclosed with the return.

If the device and/or accessories are returned to the manufacturer, the complete return shipment must be cleaned and decontaminated by the sender. If returns are not cleaned and/or decontaminated or are insufficiently cleaned and/or decontaminated, this will be performed by the manufacturer and charged to the sender.

The original transport locks must be attached for return shipment, see → *Chapter 4 'Transport and storage' on page 20*. The device must be shipped in its original packaging.

4 Transport and storage

4.1 Transport and storage conditions

Transport conditions



NOTICE

Damage to the device due to failure to use the transport locks.

- Secure the transport locks before transporting the device.



NOTICE

Damage to the device due to condensation.

There is a risk of condensation forming on electrical components when component surfaces are cold and the surrounding air is warmer. The condensation that forms may cause a short circuit and/or destroy electronics.

- Warm the device up for at least 3 hours in a warm room before connecting it to the mains.
or
- Warm up for 30 minutes in a cold room.

- Before transporting, fasten the transport lock and disconnect the device from the mains socket.
- The transport temperature must be between -20 °C and +60 °C.
- Humidity must not be condensing. Humidity must be between 10% and 80%.
- Be aware of the weight of the device.
- When transporting using a transport aid (e.g., a pallet truck), the transport aid must be able to carry at least 1.6 times the transport weight of the device.
- Secure the device to prevent it tipping over and falling down during transport.
- Never transport the device sideways or upside down.

Storage conditions

- The device must be stored in the original packaging.
- Only store the device in dry rooms.
- The storage temperature must be between -20 °C and +60 °C.
- Humidity must not be condensing. Humidity must be between 10% and 80%.

5 Commissioning

5.1 Unpacking the centrifuge

**CAUTION**

Danger of crushing due to parts falling out of the transport packaging.

- Keep the device balanced during the unpacking process.
- Only open the packaging at the points provided for this purpose.

**CAUTION**

Risk of injury from lifting heavy loads.

- Provide an adequate number of helpers.
- Note the weight. See *Chapter 3 'Device overview' on page 10*.

**NOTICE**

Damage to the device due to improper lifting.

- Do not lift the centrifuge by the control panel or the control panel holder.

Personnel:

- Trained user




Fig. 27: Unpacking

- 1 Hexagon nut
- 2 Flat
- 3 Device feet
- 4 Castors
- 5 Metal rail
- 6 Wooden beams

1. → Remove the packaging.
2. → Remove the wooden beam (6).
3. → Attach the metal rails (5) to the wooden pallet using two nails each.
4. → Slide the wooden beam (6) under the metal rails (5) to support them.
5. → Place an open-end spanner (size 10 mm) on the surfaces (2) and turn the device feet (3) upwards as far as possible.

6. Carefully roll the centrifuge off the wooden pallet over the metal rails (5).
7. Push the centrifuge to its installation location.
8. Place the open-end spanner (size 10 mm) on the flats (2) and turn the device feet (3) down until the castors (4) are no longer in contact with the ground.
9. Align the centrifuge so it is horizontal by turning the device feet (3).
10. Turn the hexagon nuts (7) upwards using the open-end spanner provided (size 19 mm) and screw them in to lock the device feet in position (3).

5.2 Setting up and connecting the centrifuge

Setting up the centrifuge



WARNING

Risk of injury due to failing to maintain a sufficient distance to the centrifuge.

- As per EN / IEC 61010-2-020, no persons, hazardous materials or objects may be present within a **safety zone of 300 mm** around the centrifuge during a centrifugation run.
- A distance of **300 mm** from the ventilation slots and ventilation openings of the centrifuge must be maintained.



CAUTION

Risk of crushing and damage to the device due to it falling down because of vibration-induced position alterations.

- Place the device on a stable and level surface.
- Select the installation surface dependent on the weight of the device.



NOTICE

Damage to the samples and the device if the ambient temperature exceeds or falls below the respective maximum/minimum permissible ambient temperature.

- Comply with the maximum and minimum permissible ambient temperatures for installation of the device.
- Do not place the device next to a heat source.
- Do not expose the device to direct sunlight.
- Do not expose the device to frost.

Personnel:

■ Trained user

1. Place the device on a stable and level surface.
2. Maintain a distance of 300 mm around the device.
3. Comply with the ambient conditions in the technical data (*→ Chapter 3.1 'Technical data' on page 10*).

Connecting the centrifuge



NOTICE

Damage to the device by unauthorised personnel

- Tampering with and modifications to devices by unauthorised persons are at the operating organisation's own risk and will result in the loss of all warranty and liability claims.



NOTICE

Damage to the device due to condensation.

There is a risk of condensation forming on electrical components when component surfaces are cold and the surrounding air is warmer. The condensation that forms may cause a short circuit and/or destroy electronics.

- Warm the device up for at least 3 hours in a warm room before connecting it to the mains.
or
- Warm up for 30 minutes in a cold room.

Personnel:

- Trained user

1. → Types 4950-08, 4950-58, 4950-78 and 4950-88 are permanently connected devices.

With permanently connected devices, a switch must be fitted in the building installation to disconnect the mains supply to the device, in accordance with the laboratory equipment standard EN / IEC 61010-1.

The switch must be located near the device, be readily accessible to the user and be marked as the means of isolation for this device.

It must be possible to secure the switch to prevent it being switched on again.

2. → A type B residual current circuit breaker must be used if the device is additionally protected with a residual current circuit breaker in the building installation.

When using a different type, the residual current circuit breaker may either not switch off the unit if there is a fault on the unit, or it may switch off the unit even though there is no fault on the unit.

3. → Types 4950-08, 4950-70, 4950-78, 4950-80 and 4950-88 must be connected in accordance with the instructions for setup and installation (AH4950).

4. → Centrifuge with PE connector:

If necessary, connect the PE connector on the back of the device to an additional medical potential equalisation system.

5. → Centrifuge with optical interface:

Connect the optical interface of the centrifuge to the PC using a fibre optic cable.

6. → Check whether the mains voltage matches the specification on the rating plate.

7. → For types 4950 and 4950-80:

Connect the device to a standard mains socket using the mains cable.

5.3 Switching the centrifuge on and off.

Switching the centrifuge on

Personnel:

- Trained user

→ Set the mains switch to */I/*.

⇒ The buttons flash, depending on the centrifuge type.

The following indicators appear one after the other, depending on the centrifuge type:

- the centrifuge model
- the last rotor code recognised by the rotor detection and the maximum rotor speed
- the program version
- When the lid is closed: ‘OPEN OEFFNEN’ indicator
- When the lid is open: The centrifugation data of the last program used or program 1.

Immediate display of centrifugation data after switching on

1. → Set the mains switch to */I/*.

2. → Press any button (except the *[STOP]* button) at the first visual change in the display (inverse display).

⇒ Centrifugation data is displayed.

Switching off the centrifuge

The rotor is stationary.

→ Set the mains switch to */O/*.

6 Operation

6.1 Opening and closing the lid

Opening the lid

Personnel:

- Trained user

The centrifuge is switched on.

The rotor is stationary.

→ Press the *[Lid]* button.

⇒ The lid unlocks by means of a motor.

The light on the *[Lid]* button goes out.

The ‘*Lid open*’ indicator appears.

Closing the lid



CAUTION

Danger of crushing when closing the lid.

Danger of crushing fingers when the closing motor pulls the lid against the seal.

- No parts of the body should be in the hazard zone of the lid when closing the lid.
- To close the lid, press on the lid from above.

**NOTICE**

Damage to the device caused by the lid slamming.

- Close the lid slowly.
- Do not slam the lid.

Personnel:

- Trained user

→ Close the lid and press down gently on the handle strip.

► The lid locks using a motor.

The *Lid*/button lights up.

The *'Lid closed'* indicator appears.

6.2 Removing and installing the rotor

Removing the rotor with a clamping nut




Fig. 28: Rotor installation and removal

- 1 Groove
- 2 Motor shaft
- 3 Driver
- 4 Hole

Personnel:

- Trained user

1. → Open the lid.

2. → Loosen the rotor clamping nut using the supplied spanner.

► After passing the working point for lifting the rotor, the rotor detaches from the cone of the motor shaft (2).

3. → Turn the clamping nut until the rotor can be lifted off the motor shaft.

4. → Remove the rotor.

Installing the rotor with a clamping nut

Personnel:

- Trained user

The lid is open.

1. → Clean the motor shaft (2) and rotor hole (4).

2. → Lightly grease the motor shaft (2), see ► Chapter 8.2 'Cleaning and disinfection instructions' on page 43.

3. → Place the rotor vertically on the motor shaft (2).

The driver (3) of the motor shaft must be in the groove (1) of the rotor. The orientation of the groove is marked on the rotor.

4. → Hand-tighten the rotor clamping nut using the supplied spanner.

5. → Check that the rotor is firmly seated.

6.3 Inserting and removing buckets

Inserting buckets



NOTICE


Damage to the device due to imbalances caused by incorrect loading of the rotor.

- Load all swing-out rotor locations with the same buckets.



Buckets marked with the number of the rotor location may only be used there.

Buckets marked with a set number may only be used together.




1. Check that the rotor is firmly seated.
2. Grease the trunnions (3).
3. Insert the bucket (1) into the rotor from above. The trunnions (3) must be in the grooves (2).
4. Push the bucket (1) down as far as it will go.

Removing the bucket

- Pull the bucket (1) vertically upwards out of the rotor.

6.4 Inserting and removing adapters

Inserting



the adapter

- Insert the adapter (1) vertically into the bucket (2) from above.

removing

- Remove the adapter (1) vertically upwards out of the bucket (2).

6.5 Loading

Filling centrifuge tubes



WARNING

Risk of injury from contaminated sample material.

Contaminated sample material escapes from the sample tube during centrifugation.

- Use centrifuge tubes with special screw caps for hazardous substances.
- For risk group 3 and 4 materials, use a biosafety system in addition to the sealable centrifuge tubes (see WHO's 'Laboratory Biosafety Manual').



NOTICE

Damage to the device due to highly corrosive substances.

Highly corrosive substances may impair the mechanical strength of rotors, buckets and accessories.

- Do not centrifuge highly corrosive substances.



Standard glass centrifuge tubes can be loaded up to RCF 4000 (DIN 58970 part 2).

Personnel:

- Trained user

→ Fill centrifuge tubes outside the centrifuge.

The maximum capacity of the centrifuge tubes specified by the manufacturer must not be exceeded.


With angle rotors, the centrifuge tubes must only be filled to the extent that no liquid can be ejected from the tubes during the centrifugation run.

It must be ensured that there is a uniform fill level in the tubes in order to keep the weight differences in the centrifuge tubes as low as possible.

Loading the swing-out rotors

Personnel:

- Trained user



The following must be observed when using blood bags:

1. Check that the rotor is firmly seated.
2. The centrifuge tubes must be distributed symmetrically across all rotor locations.

The weight of the permissible filling capacity is indicated on each rotor. The weight must not be exceeded.

No liquid must be allowed to enter the buckets and the centrifuging chamber when loading the buckets and swinging them out during the centrifugation run.

For containers with rubber inserts, there must always be the same number of rubber inserts under the centrifuge tubes.

All rotor locations must be filled with the same buckets. Certain buckets are marked with the number of the rotor location. The buckets must only be inserted in the corresponding rotor location.


Buckets marked with a set number (for example S001/4) must only be used in the set.

1. The differences can be compensated with balancing weights if the buckets are not filled with the same weight.
2. Empty buckets can be loaded with compensating inserts if there are insufficient blood bag systems available to load the rotor fully.
3. If necessary, fine-tune with the tare weights supplied.

Loading the angle rotors

Personnel:

- Trained user






1. → Check that the rotor is firmly seated.
2. → The centrifuge tubes must be distributed evenly over all locations on the rotor.

No liquid must be allowed to enter the rotor and the centrifuging chamber when loading the rotor.

With rotors, the centrifuge tubes must only be filled to the extent that no liquid can be ejected from the tubes during the centrifugation run.

The weight of the permissible filling capacity is indicated on each rotor. The weight must not be exceeded.



6.6 Opening and closing the biosafety system

6.6.1 Explanation

The user must take appropriate actions when centrifuging hazardous substances or mixtures of substances that are toxic, radioactive or contaminated with pathogenic microorganisms.

Centrifuge tubes with special screw caps for hazardous substances must always be used.

For materials of risk group 3 and 4, a biosafety system must be used in addition to the sealable centrifuge tubes (see the World Health Organisation's "Laboratory Biosafety Manual").

In a biosafety system, a bioseal (sealing ring) prevents droplets and aerosols from escaping.

If the bucket of a biosafety system is used without the lid, the sealing ring must be removed from the bucket to prevent damage to the sealing ring during the centrifugation run.

Damaged biosafety systems are no longer microbiologically tight.

If no biosafety system is used, a centrifuge is not microbiologically tight in the sense of the EN / IEC 61010-2-020 standard.

Storage of biosafety systems

Biosafety systems must only be stored with the lid open to avoid damage to the sealing rings during storage.

6.6.2 Lid with screw cap and hole




Fig. 29: Biosafety system

- 1 Rotary handle
- 2 Rotor
- 3 Lid
- 4 Key

Closing

1. Place the lid (3) centrally on the rotor (2).
2. Insert the supplied key (4) into the hole in the rotary handle (1).
3. Turn the lid (3) at the key (4) clockwise until it is tightly closed.

Opening

1. Insert the supplied key (4) into the hole in the rotary handle (1).
2. Turn the lid (3) at the key (4) anticlockwise until it is open.
3. Remove the lid (3) from the rotor (2).

6.6.3 Lid with screw cap




Fig. 30: Biosafety system

- 1 Lid
- 2 Bucket

Closing

1. Place the lid (1) centrally on the bucket (2).
2. Turn the lid (1) clockwise until it is tightly closed.

Opening


1. Turn the lid (1) anticlockwise until it is open.
2. Remove the lid (1) from the bucket (2).

6.7 Packing instructions, HettLiner

Packing before centrifugation



Ensure that the plastic insert cannot tip over when loading and unloading the inserts (if necessary, use loading aid 4509).



1. → Insert the blood bag (1) into the insert (2).



2. → Hold the blood bag by the connectors (1) and push the support plate (2) on the outer side of the blood bag into the insert from top to bottom.

Ensure that the lower edge of the support plate stands on the floor as completely as possible.



3. → Fold the support plate outwards and press it down until the folded edge of the support plate is level with the liquid level of the blood bag.

The upper edge of the support plate must not project too far out of the insert during centrifugation due to the risk of its becoming caught in the rotor arms.

Note the position of the loop (1) so that it can be reached after centrifugation.

- 4.** If available, fold empty satellite bag(s) and pack differently depending on the appropriate accessories and filling volume of the blood bag. It is advantageous to fold the satellite bags and pack them on the outside between the folded-down support plate and the outer wall of the insert.
Ensure that the silicone plate does not slip when doing this.
If necessary, the silicone plate can be retained and held against the loop when packing the satellite bag.
The position of the loop must be checked afterwards.
- 5.** Place the connections over the support plate so that the valves cannot break.
Ensure that the hoses do not project out of the insert.
Stow hose sections projecting over the edge of the insert between the folded-down support plate and the insert wall.
- 6.** If necessary, balance weights should be placed between the folded-down support plate and the receptacle wall.

Unpacking after centrifugation

- 1.** Pull the satellite bag out of the insert while holding the silicone plate in place with one hand.
- 2.** Slowly pull out the folded-down part of the support plate at the loop provided for this purpose.
Return the support plate to its original shape in a controlled manner. The folded-down part of the support plate may spring back and mix blood components.
- 3.** Remove the remaining blood bag either together with the support plate or after removing the support plate from the insert.

6.8 Centrifugation

6.8.1 Centrifugation in continuous operation

Personnel:

■ Trained user

- 1.** Press the *[/t]* button repeatedly until the input field of the '*t/min:*' parameter is reverse highlighted.
- 2.** Select the value 0 using the rotary knob.
- 3.** Press the *[/t]* button repeatedly until the input field of the '*/:sec*' parameter is reverse highlighted.
- 4.** Select the value 0 using the rotary knob.
→ The input field displays '*---:--*'.
- 5.** Press the *[START]* button.
→ The centrifugation run is started.
The '*Rotation*' indicator is lit up while the rotor is turning.
The timing starts at 00:00.
The rotor speed or the resulting RCF value, the temperature in the centrifuging chamber and the elapsed time are displayed during the centrifugation run.
- 6.** Press the *[STOP]* button to cancel the centrifugation run.
Ramp-down takes place with the set ramp-down parameters.
→ '*OPEN OEFFNEN*' is displayed.

6.8.2 Centrifugation with time preselection

Personnel:

- Trained user
1. → Press the *[f]* button repeatedly until the input field of the 't/min:' parameter is reverse highlighted.
 2. → Use the *[Rotary knob]* to set the desired value.
 3. → Press the *[f]* button repeatedly until the input field of the 't/:sec' parameter is reverse highlighted.
 4. → Use the *[Rotary knob]* to set the desired value.
 5. → Press the *[START]* button.
 - The centrifugation run is started.
- The '*Rotation*' indicator is lit up while the rotor is turning.
- The rotor speed or the resulting RCF value, the temperature in the centrifuging chamber and the remaining time are displayed during the centrifugation run.
6. → Ramp-down takes place with the selected ramp-down parameters after the time has elapsed or if the centrifugation run is cancelled by pressing the *[STOP]* button.
 - '*OPEN OEFFNEN*' is displayed.

6.8.3 Changing settings during centrifugation

The runtime, speed, relative centrifugal force (RCF/RZB), ramp-up and ramp-down parameters, and the temperature (only for devices with cooling) can be changed during centrifugation.

These parameters can only be changed individually and one after the other.

1. → Change the value of the desired parameter using the *[Rotary knob]*
 2. → Press the *[START]* button.
 - The values of the current program are copied to program location '----' and updated with the changed value.
- The original program is not overwritten.

6.9 Quick stop function

Personnel:

- Trained user
- Press the *[STOP]* button twice.
- The '*STOP*' indicator flashes.

Ramp-down with brake level "R9" (shortest ramp-down time) is displayed and executed.

If brake level "R0" was selected, the ramp-down time is longer than with brake level "R9" for technical reasons.

7 Software operation

7.1 Key switches

The keys must be stored in such a way that they are protected from unauthorised access.

Key position	Function
Left key position	'LOCK 1' is displayed. Programs can only be retrieved, not changed.
Right key position	'LOCK 2' is displayed. No programs can be retrieved or changed.
Middle key position	no status indicator. No program lock. Programs can be retrieved and changed.

7.2 Centrifugation parameters

7.2.1 Ramp-up and ramp-down parameters

The set ramp-up and ramp-down parameters are displayed.



x: 1-9 = ramp-up level, t = ramp-up time

y: R1-R9, B1-B9 = brake level, R0 = unbraked ramp-down, t = ramp-down time, n(*) = brake cut-off speed

Ramp-up level

1. Press the [Ramp-up parameters] button repeatedly until the 'Ramp-up level' parameter or the 'Ramp-up time' parameter is displayed and the input field is reverse highlighted.

2. Use the [Rotary knob] to set the desired level.

Ramp-up time

1. Press the [Ramp-up parameters] button repeatedly until the 'Ramp-up time min:sec' parameter is displayed and the input field is reverse highlighted.

2. Use the [Rotary knob] to set the desired level.

If a ramp-up time is set that is longer than the runtime, the centrifugation run ends before the set speed is reached.

Brake level

1. Press the [Ramp-down parameters] button repeatedly until the 'Ramp-down level' parameter or the 'Ramp-down time' parameter is displayed and the input field is reverse highlighted.

2. Use the [Rotary knob] to set the desired level.

B-brake levels can only be set for special rotors.

Ramp-down time

No ramp-down time can be set if a brake cut-off speed is set.

1. Press the [Ramp-down parameters] button repeatedly until the 'Ramp-down level min:sec' parameter is displayed and the input field is reverse highlighted.

2. Use the [Rotary knob] to set the desired level.

Brake cut-off speed

1. Press the [Ramp-down parameters] button repeatedly until the 'n*/RPM' parameter is displayed and the input field is reverse highlighted.

2. Use the [Rotary knob] to set the desired level.

7.2.2 Runtime



For continuous operation, the minutes and seconds must be set to zero.

Continuous operation is indicated in the indicator by means of the '---:--' symbol.

1. → Press the *[f]* button repeatedly until the input field of the 't/min:' parameter is reverse highlighted.
2. → Use the *[Rotary knob]* to set the desired value.
3. → Press the *[f]* button repeatedly until the input field of the 't/sec' parameter is reverse highlighted.
4. → Use the *[Rotary knob]* to set the desired value.

7.2.3 Speed, RPM

1. → Press the *[n]* button repeatedly until the 'RPM' parameter is displayed and the input field is reverse highlighted.
2. → Use the *[Rotary knob]* to set the desired value.

Maximum rotor speed indicator

1. → Press the *[n]* button repeatedly until the 'RPM' parameter is displayed and the input field is reverse highlighted.
2. → Press and hold the *[n]* button.
 - The maximum rotor speed (n-max rotor) is displayed.

7.2.4 Integral RCF

Integral RCF is a measure of the sedimentation effect ($\int n^2 dt$). The value is used to compare centrifugation runs.

- Press and hold the *[Integral RCF]* button.
 - 'Integral RCF' is displayed.

7.2.5 temperature

1. → Press the *[Temperature and centrifuging radius]* button repeatedly until the 'T/C°' parameter is displayed and the input field is reverse highlighted.
2. → Use the *[Rotary knob]* to set the desired value.

7.2.6 Relative centrifugal force, RCF

The relative centrifugal force RCF is dependent on the speed and the centrifuging radius.

The relative centrifugal force RCF is stated as a multiple of the acceleration due to gravity (g).

The relative centrifugal force RCF is a dimensionless numerical value and is used to compare the separation and sedimentation performance.

$$RCF = \left(\frac{RPM}{1000} \right)^2 * r * 1,118$$

$$RPM = \sqrt{\frac{RCF}{r * 1,118}} * 1000$$

RCF = Relative Centrifugal Force

RPM = speed

r = centrifuging radius in mm = distance from the centre of the axis of rotation to the bottom of the centrifuge tube.

7.2.7 Setting the relative centrifugal force (RCF/RZB)

1. → Press the *[RCF]* button repeatedly until the 'RCF/RZB' parameter is displayed and the input field is reverse highlighted.
2. → Use the *[Rotary knob]* to set the desired value.

Maximum rotor RCF indicator

1. → Press the *[RCF]* button repeatedly until the 'RCF/RZB' parameter is displayed and the input field is reverse highlighted.
2. → Press and hold the *[RCF]* button.
→ The maximum RCF of the rotor (RCF-max rotor) is displayed.

7.2.8 Centrifugation of substances or mixtures of substances with a density higher than 1.2 kg/dm³

The density of the substances or mixtures of substances must not exceed 1.2 kg/dm³ during centrifugation at maximum speed. The speed must be reduced for substances or substance mixtures with a higher density. The permissible speed can be calculated using the following formula:

$$\text{Reduzierte Drehzahl } (n_{red}) = \sqrt{\frac{1,2}{\text{höhere Dichte (kg/dm}^3)}} * \text{maximale Drehzahl (RPM)}$$

For example: Maximum speed 4000 RPM, density 1.6 kg/dm³

$$n_{red} = \sqrt{\frac{1,2(\text{kg/dm}^3)}{1,6(\text{kg/dm}^3)}} * 4000 \text{ RPM} = 3464 \text{ RPM}$$

If, in exceptional cases, the maximum load indicated on the bucket is exceeded, the speed must also be reduced. The permissible speed can be calculated using the following formula:

$$\text{Reduzierte Drehzahl } (n_{red}) = \sqrt{\frac{\text{maximale Beladung (g)}}{\text{tatsächliche Beladung (g)}}} * \text{maximale Drehzahl (RPM)}$$

For example: Maximum speed 4000 RPM, maximum load 300 g, actual load 350 g

$$n_{red} = \sqrt{\frac{300 \text{ g}}{350 \text{ g}}} * 4000 \text{ RPM} = 3703 \text{ RPM}$$

Please contact the manufacturer if you are not sure.

7.2.9 Centrifuging radius

1. → Press the *[Temperature and centrifuging radius]* button repeatedly until the 'r/mm' parameter is displayed and the input field is reverse highlighted.
 2. → Use the *[Rotary knob]* to set the desired value.
- Changing the radius, automatically adjusts the RCF/RZB value: this is indicated by flashing.

7.3 Programming

7.3.1 Opening or loading programs

1. → Use the *[PROG]* button to select the 'PROG no.' parameter. The input field is reverse highlighted.
2. → Use the *[Rotary knob]* to set the desired program location.

3. ➤ Press the *[RCL]* button.
 - The centrifugation data of the desired program location is displayed.

7.3.2 Entering or changing programs

1. ➤ Set the desired parameters.
2. ➤ Use the *[PROG]* button to select the 'PROG no.' parameter. The input field is reverse highlighted.
3. ➤ Use the *[Rotary knob]* to set the desired program location.
If the program location indicator flashes, this program location is already assigned centrifugation data. In this case, set a free program location, or overwrite the centrifugation data by continuing.
4. ➤ Press the *[STO]* button.
 - Settings are stored in the desired program location.
5. ➤ Press the *[STO]* button twice.
 - Centrifugation data that is already stored will be overwritten.

7.3.3 Automatic buffer

The buffer includes the program locations "----" and 90 to 99. Changed centrifugation data is automatically saved to program location "----" every time a centrifugation run is started. The changed centrifugation data, of the last 11 centrifugation runs, are stored in the buffer and can be retrieved.

7.4 Rotor detection

- Rotor detection is performed after starting a centrifugation run.
- If the rotor has been changed, the centrifugation run is cancelled after rotor detection. The rotor code (R) and the maximum rotor speed (n-max) of the newly detected rotor are displayed.
- If the maximum speed of the rotor used is less than the set speed, the speed is limited to the maximum rotor speed.

7.5 Cooling (for centrifuges with cooling)

7.5.1 Instructions, cooling

For centrifuges with the heating/cooling option, the temperature setpoint can be adjusted from -20 °C to +60 °C. If the actual temperature deviates from the set temperature by more than 5 °C, this is signalled by a flashing temperature value indicator.

The lowest achievable temperature is rotor dependent.

7.5.2 Standby cooling

The centrifuging chamber is cooled to the preselected temperature when the rotor is at a standstill and the lid is closed. The display shows the temperature setpoint.

7.5.3 Precooling the rotor

For rapid precooling of the unloaded rotor and accessories, we recommend a centrifugation run with the continuous operation settings and a speed of approx. 20% of the maximum rotor speed.

7.6 Heating (for centrifuges with heating)

During the centrifugation run, the centrifuging chamber is heated to the preselected temperature if required. The heating is switched off when the rotor is at a standstill.



CAUTION

Danger of burns from hot surfaces.

The surface temperature of the heating element in the centrifuging chamber can be up to 500 °C or 932 °F.

- Do not touch the heating element.



NOTICE

Damage to plastic buckets due to excessive temperature

- Plastic buckets may only be used at temperatures up to a maximum of 40 °C or 104 °F.

7.7 Machine Menu

7.7.1 Querying system information

The following system information can be queried:

- Centrifuge model
- Maximum speeds of the various rotor codes
- Centrifuge program version
- Frequency converter type
- Program version for the frequency inverter

The rotor is stationary.

1. Press and hold the *[/t]* button.
→ The audible signal ‘SOUND / BELL’ [sounds] after 8 seconds.
2. Press the *[/t]* button.
→ The operating hours ‘CONTROL.’ are displayed.
3. Press the *[/t]* button.
→ The date and time are displayed.
4. Press the *[/t]* button.
→ The machine and cooling version ‘VERS 12 °C /* 03’ is displayed.
5. Press the *[/t]* button.
→ The operating hours of the frequency converter ‘FC/CCI XX h’ are displayed.
6. Press the *[/t]* button.
→ The type of frequency converter ‘FU/CCI’ is displayed.
7. Press the *[/t]* button.
→ The program version of the frequency converter ‘FU/CCI – S.’ is displayed.
8. Press the *[/t]* button.
→ The program version of the power supply circuit board ‘°C /* – S. 01.07’ is displayed.
9. Press the *[STOP/OPEN]* button to exit the menu

7.7.2 Querying operating hours

The rotor is stationary.

1. → Open the lid.
2. → Press and hold the *[/t]* button.
 - ‘SOUND / BELL XXX’ is displayed after 8 seconds.
3. → Press the *[/t]* button.
 - ‘CONTROL:’ and the operating hours are displayed.

The operating hours indicator goes off automatically after 10 seconds.

7.7.3 Audible signal

7.7.3.1 General

The audible signal sounds after the following settings:

OFF	<ul style="list-style-type: none">■ after a problem occurs in the 2 s interval.
ON1	<ul style="list-style-type: none">■ after a problem occurs in the 2 s interval.■ after completion of the centrifugation run and rotor standstill in the 30 s interval.
ON2	<ul style="list-style-type: none">■ after a problem occurs in the 2 s interval.■ after completion of the centrifugation run and rotor standstill in the 30 s interval.■ every time the button is pressed.

Opening the lid or pressing any button stops the audible signal.

7.7.3.2 Setting an audible signal

1. → Open the lid.
2. → Press and hold the *[/t]* button.
 - ‘SOUND / BELL ON1’, ‘SOUND / BELL ON2’ or ‘SOUND / BELL OFF’ is displayed after 8 seconds.
3. → Use *[Rotary knob]* to set ‘OFF’, ‘ON1’ or ‘ON2’.
4. → Press the *[START]* button.
 - The setting is stored.

‘*** OK ***’ is displayed briefly.

7.7.4 Centrifugation data displayed after switching on

The centrifugation data of program 1 or the last program used is displayed after switching on.

1. → Set the mains switch to *///*.
2. → Press the *[STOP]* button at the first visual change in the display (inverse display).
 - ‘PROGRAM 1, LAST PROGRAM’ is displayed.
3. → Use the *[Rotary knob]* to set the desired function.

4. Press the *[START]* button.
⇒ The settings are stored.
‘***OK***’ is displayed briefly.

7.7.5 Setting the date and time

- The rotor is stationary.
1. Open the lid.
 2. Press and hold the *[t]* button.
⇒ ‘SOUND / BELL’ is displayed after 8 seconds.
 3. Press the *[t]* button twice.
⇒ The date and time are displayed
a: Year
mon: Month
d: Day
h: Hours
min: Minutes
 4. Press the *[Temperature and centrifuging radius]* button repeatedly until the desired parameter is displayed and the input field is reverse highlighted.
 5. Use the *[Rotary knob]* to set the desired value.
 6. Press the ‘Start’ button.
⇒ The settings are stored.
‘***OK***’ is displayed briefly.

7.8 Program links

7.8.1 Linking programs or changing a program link



Program linking is only possible with programs where ramp-up and brake levels are set.

The programs must be saved in the desired order before linking, by means of either program entry or program retrieval.

The program locations must be consecutive (e.g. program locations 10+11+12).


Linking programs

1. Use the *[PROG]* button to select the ‘PROG no.’ parameter. The input field is reverse highlighted.
2. Use the *[Rotary knob]* to set the program location of the initial program (XX+).
3. Press the *[RCL]* button.
⇒ The centrifugation data of the desired program location is displayed
4. Press the *[PROG]* button twice.
⇒ Parameter PR-PART is selected.
The input field is reverse highlighted.

5. → Press the *[STO]* button twice.
 - The program is linked and the program number of the next program location (+XX+) is displayed.
 6. → Press the *[RCL]* button twice.
 - The centrifugation data of the desired program location is displayed
 7. → Press the *[STO]* button twice.
 - The program is linked and the program number of the next program location (+XX+) is displayed.
 8. → Repeat the last two steps until all programs are linked.
 9. → Press the *[PROG]* button.
 - The program number of the end program (+XX) is displayed.
- Changing a program link**
1. → Open the desired program.
 2. → Change the desired parameter.
 3. → Save the changed centrifugation data to the same program location again.
 - Saving removes the program link.
 4. → Link programs again.

7.8.2 Centrifugation run with program link

1. → Press the *[PROG]* button twice.
 - Parameter PR-PART is selected.
The input field is reverse highlighted.
2. → Use the *[Rotary knob]* to set the program location of the initial program (XX+).
3. → Press the *[RCL]* button.
 - The centrifugation data of the desired program location is displayed
4. → Press the *[START]* button.
 - The centrifugation run is started.
The ‘Rotation’ indicator appears as long as the rotor is turning.
The ramp-up and brake levels of the program link are displayed.
 - Initial program (XX+)
x: Ramp-up level of the initial program x
 - Follow-up program (+XX+)
x: Ramp-up level of the follow-up program x
 - End program (+XX)
x: Ramp-up level of the end program
y: Brake level of the end program
5. → Ramp-down takes place with the brake level of the end program after the time in the end program has elapsed.
If the centrifugation run is cancelled by pressing the *[STOP]* button, ramp-down takes place at the brake level of the program currently running.



7.8.3 Deleting program links

- 1.** Use the *[PROG]* button to select the 'PROG no.' parameter. The input field is reverse highlighted.
- 2.** Use the *[Rotary knob]* to set the program location of the initial program (XX+).
- 3.** Press the *[RCL]* button.
 - ⇒ The centrifugation data of the desired program location is displayed
- 4.** Press the *[PROG]* button twice.
 - ⇒ Parameter 'PR-PART' is displayed.
The input field is reverse highlighted.
- 5.** Press the *[STO]* button twice.
- 6.** Press the *[PROG]* button.


8 Cleaning and care

8.1 Overview table

Chap.	Task to execute	if required	daily	weekly	Annually	Page
8	Cleaning and care					42
8.3	Cleaning					43
8.3	Cleaning the device		X			43
8.3	Cleaning the biosafety systems			X		44
8.3	Cleaning the accessories			X		44
8.4	Disinfection					44
8.4	Disinfecting the device	X				44
8.4	Disinfecting the accessories	X				44
8.5	Maintenance					45
8.5	Greasing the rubber seal of the centrifuging chamber			X		45
8.5	Greasing the rubber seal of the biosafety system			X		45
8.5	Trunnion greasing			X		45
8.5	Checking the accessories			X		45
8.5	Checking the biosafety system			X		45
8.5	Inspecting the centrifuging chamber for damage				X	45
8.5	Greasing the motor shaft				X	45

Chap.	Task to execute	if required	daily	weekly	Annually	Page
8.5	Accessories with a limited service life	X				46
8.5	Calculation of the number of run cycles performed	X				46
8.5	Replacing centrifuge tubes	X				46

8.2 Cleaning and disinfection instructions



- The device and its accessories must not be cleaned in dishwashers.
- Only perform hand cleaning and liquid disinfection.
- The water temperature must not exceed 25 °C.
- To prevent any corrosion due to use of detergents or disinfectants, it is essential to follow the special application instructions provided by the manufacturers of the detergent or disinfectant.

Disinfectant:

- Surface disinfectant (not disinfectant for hands or instruments)
- Ethanol as the sole active substance.
Do not use an ethanol-propanol mixture to disinfect the viewing window in the lid of the device.
- Concentration is not less than 30 %
- pH: 6 – 8
- Non-corrosive

8.3 Cleaning

Cleaning the device

1. → Open the lid.
2. → Switch off the device and disconnect it from the power supply.
3. → Remove accessories.
4. → Clean the centrifuge housing and the centrifuging chamber with soap or a mild detergent and a damp cloth.
5. → Remove any detergent residues with a damp cloth after using detergents.
6. → The surfaces must be dried immediately after cleaning.
7. → Dry the centrifuging chamber with an absorbent cloth if condensation forms.

Cleaning the biosafety systems

- 1.** Clean the biosafety system using the detergent and a damp cloth.
- 2.** Remove any detergent residues with a damp cloth after using detergents.
- 3.** Dry the accessories immediately after cleaning using a lint-free cloth and oil-free compressed air. Dry all cavities completely using oil-free compressed air.

Cleaning the accessories

- 1.** Clean the accessories using the detergent and a damp cloth.
- 2.** Remove any detergent residues with a damp cloth after using detergents.
- 3.** Dry the accessories immediately after cleaning using a lint-free cloth and oil-free compressed air. Dry all cavities completely using oil-free compressed air.

8.4 Disinfection



Disinfection must always be preceded by cleaning of the components concerned.

See ➔ Chapter 8 'Cleaning and care' on page 42



Disinfectant concentration and application time according to the manufacturer's instructions.

Disinfecting the device



CAUTION

Risk of injury due to ingress of water or other liquids.

- Protect the device against external liquids.
- Do not disinfect the device using spray.

- 1.** Open the lid.
- 2.** Switch off the device and disconnect it from the power supply.
- 3.** Remove accessories.
- 4.** Clean the housing and centrifuging chamber using disinfectant.
- 5.** Remove any disinfectant residues with a damp cloth after using disinfectants.
- 6.** The surfaces must be dried immediately after cleaning.

Disinfecting the accessories

- 1.** Disinfect the accessories using the disinfectant.
- 2.** Wet all cavities with bubble-free disinfectant.
- 3.** Remove the disinfectant residues or leave them to dry after using disinfectants.

Autoclaving

The following accessories may be autoclaved at 121 °C / 250 °F (20 min):

- Swing-out rotors
- Aluminium angle rotors
- Metal buckets

- Lid with bioseal
- Inserting

No statement can be made about the resulting degree of sterility.

The lids of the rotors and bucket must be removed before autoclaving.

Autoclaving accelerates the ageing of materials. It may cause changes to colours. After autoclaving, the rotors and accessories are to be visually inspected for damage and any damaged parts are to be replaced immediately.

The sealing ring in question is to be replaced if there are signs of cracking, embrittlement or wear. For lids with non-replaceable sealing rings, the whole lid must be replaced.

The sealing rings must be replaced after autoclaving to ensure the tightness of the biosafety systems.

8.5 Maintenance

Greasing the rubber seal of the centrifuging chamber

→ Rub the sealing ring lightly with a rubber care product.

Greasing the rubber seal of the biosafety system

→ Rub the sealing ring lightly with a rubber care product.

Trunnion greasing

1. → Remove accessories.
2. → Clean the trunnions.
3. → Remove any detergent residues with a damp cloth after using detergents.
4. → Grease the trunnions and suspension with Hettich Tubenfett 4051.
5. → Excess grease in the centrifuging chamber must be removed.

Checking the accessories

1. → The accessories are to be checked for wear and corrosion damage.
2. → Check that the rotor is firmly seated.

Checking the biosafety system

1. → Visually check all parts of the biosafety system for damage.
2. → Check the correct installation position of the sealing ring(s) of the biosafety system.
3. → Replace the damaged parts of the biosafety system.
4. → Replace the sealing ring in question immediately if there are signs of cracking, embrittlement or wear. For lids with non-replaceable sealing rings, the whole lid must be replaced.

Inspecting the centrifuging chamber for damage

→ Check the centrifuging chamber for damage.

Greasing the motor shaft

1. → Remove accessories.
2. → Clean the motor shaft.
3. → Remove any detergent residues with a damp cloth after using detergents.
4. → Grease the motor shaft with Hettich Tubenfett 4051.
5. → Excess grease in the centrifuging chamber must be removed.

Accessories with a limited service life

The use of certain accessories is time-limited. For safety reasons, the accessories must no longer be used when either the maximum number of permissible run cycles marked on them or the expiry date marked on them has been reached.

- The maximum permissible number of run cycles or the expiry date can be seen marked on the accessories.
- The centrifuge is equipped with a cycle counter.

Calculation of the number of run cycles performed

The runtime per centrifugation run and the operating hours of the device are required to calculate the number of run cycles performed (centrifugation runs). Querying the operating hours, see **Chapter 7.7 'Machine Menu' on page 38**.

The shortest runtime is to be used for the calculation if centrifugation runs with different runtimes have been carried out.

The number of run cycles performed (centrifugation runs) is calculated as follows:

Number of run cycles performed = operating hours [h] x 60 / runtime [min]

e.g.: Operating hours 2000 h, runtime 5 min

Number of run cycles performed = $2000 \times 60 / 5 = 24000$

Replacing centrifuge tubes



CAUTION

Risk of injury from broken glass.

Broken glass may cause glass splinters and contaminated liquids to be found inside the centrifuge.

- Wear cut-resistant gloves.
- Wear protective goggles and a face mask.

Broken parts of the tube, glass splinters and spilled centrifuge material must be removed completely in the event of leakage or if a centrifuge tube breaks. Glass splinters that are not removed will cause further glass breakage.

The rubber inserts and the plastic sleeves of the rotors must be replaced after a glass breakage.

Disinfection must be carried out if the material is infectious.

9 Troubleshooting

9.1 Fault description

Customer service must be notified if the fault cannot be rectified based on the fault table. State the centrifuge type and serial number. Both numbers can be seen on the type plate of the centrifuge.

* Error number does not appear on the display.

Fault description	Cause	Remedy
no display	No power. Triggering of the overcurrent protection fuse.	<ul style="list-style-type: none"> ■ Check the supply voltage. ■ The mains switch is in switch position <i>/ /</i>
TACHO - ERROR 01, 02	Tacho defective. Motor, inverter, electronics defective.	<ul style="list-style-type: none"> ■ Open the lid. ■ Set the mains switch to <i>/O/</i>.

Fault description	Cause	Remedy
TACHO - ERROR 01, 02	Tacho defective. Motor, inverter, electronics defective.	<ul style="list-style-type: none"> ■ Wait at least 10 seconds. ■ Turn the rotor vigorously by hand. ■ Set the mains switch to <i>[0]</i>. The rotor must rotate while switching on.
IMBALANCE / UNWUCHT	The rotor is unevenly loaded.	<ul style="list-style-type: none"> ■ Open the lid or hatch. ■ Check the loading of the rotor. ■ Repeat the centrifugation run.
CONTROL - ERROR 04, 06-09	Lid lock error.	<ul style="list-style-type: none"> ■ Perform a MAINS RESET.
N > MAX 05	Overspeed error	<ul style="list-style-type: none"> ■ Perform a MAINS RESET.
N < MIN 13	Underspeed error.	<ul style="list-style-type: none"> ■ Perform a MAINS RESET.
ROTORCODE 10	Rotor coding error.	<ul style="list-style-type: none"> ■ Perform a MAINS RESET.
MAINS INTERRUPT	Loss of mains power during the centrifugation run. The centrifugation run was not completed.	<ul style="list-style-type: none"> ■ Open the lid. ■ Press the <i>[START]</i> button. ■ If required: Repeat the centrifugation run.
VERSION-ERROR 12	No conformity of the electronic components, error/defect in electronics.	<ul style="list-style-type: none"> ■ Perform a MAINS RESET.
SER I/O - ERROR 30-38	Error/defect in interface.	<ul style="list-style-type: none"> ■ Perform a MAINS RESET.
° C * - ERROR 50-56, 58	Error/defect in cooling.	<ul style="list-style-type: none"> ■ Perform a MAINS RESET.
LOCK - ERROR 57	Error/defect in program lock.	<ul style="list-style-type: none"> ■ Perform a MAINS RESET.
FU/CCI-ERROR 60-83	Error/defect in motor control.	<ul style="list-style-type: none"> ■ Perform a MAINS RESET.
CONTROL - ERROR 26, 90-95, 97 - 99	Error/defect in control section.	<ul style="list-style-type: none"> ■ Perform a MAINS RESET.
N > ROTOR MAX 96	<p>Speed in the selected program greater than the maximum rotor speed.</p> <p>The rotor has been changed. The built-in rotor has a higher maximum speed than the previously used rotor. The rotor has not yet been recognised by the rotor detection.</p>	<ul style="list-style-type: none"> ■ Check and correct the speed. ■ Set a speed up to the maximum speed of the previously used rotor. Press the <i>[START]</i> button to perform rotor detection.
The entire display lights up.	-	<ul style="list-style-type: none"> ■ Notify customer service.

9.2 Perform a MAINS RESET

1. Set the mains switch to *[0]*.
2. Wait 10 seconds.

3. Set the mains switch to */I*.

9.3 Emergency release

The lid cannot be unlocked by the motor in the event of a power failure. Emergency unlocking by hand must be performed.



WARNING

Risk of electric shock due to maintenance and servicing work on live device.

- Disconnect the device from the mains before carrying out repairs and maintenance.



WARNING

Danger of cutting and crushing due to moving rotor.

- Do not open the lid until the rotor has stopped.

Personnel:

- Trained user

1. Look through the window in the lid to ensure that the rotor is stationary.
2. Insert the hex key horizontally into the hole (1) and turn clockwise until the lid opens.
3. Remove the hex key from the hole (1).
4. When the power is restored, press the *[Lid]* button so that the motorised cover lock returns to the home position (open).




Fig. 31: Emergency release

1 Hole

10 Disposal

10.1 General instructions



The device can be disposed of via the manufacturer.

A Return Material Authorisation (RMA) form must always be requested for a return.

If necessary, contact the Technical Service Department of the manufacturer.

- Andreas Hettich GmbH & Co. KG
- Föhrenstrasse 12
- 78532 Tuttlingen, Germany
- Phone: +49 7461 705 1400
- E-mail: service@hettichlab.com

**WARNING**

Risk of pollution and contamination for people and the environment.

When disposing of the centrifuge, people and the environment may be polluted or contaminated by incorrect or improper disposal.

- Removal and disposal may be carried out only by a trained and authorised service personnel.

The device is intended for the commercial sector ("Business to Business" - B2B).

According to Directive 2012/19/EU, the devices may no longer be disposed of with household waste.

The devices are assigned to the following groups according to the Stiftung Elektro-Altgeräte Register (EAR (German foundation under civil law)):

- Group 1 (heat exchangers)

The crossed-out wheelie bin symbol indicates that the device must not be disposed of with household waste. Regulations governing disposal of such devices may differ in individual countries. If necessary, contact the supplier.



Fig. 32: Household waste ban

11 Index

A	
Akustisches Signal	
aktivieren/deaktivieren.	39
Allgemeine Sicherheitshinweise.	8
An- und Auslaufparameter.	34
Anschließen der Zentrifuge.	23
Aufstellen der Zentrifuge.	22
Auspicken.	21
Ausschalten.	24
Autoklavieren.	44
B	
Befüllen.	27
Beladen.	27
Betriebsstunden	
abfragen.	39
Bio-Sicherheitssystem	
prüfen.	45
reinigen.	44
D	
Datums- und Uhreinstellung.	40
Dauerlauf.	32
Desinfektion.	44
Drehzahl RPM.	35
E	
Einschalten.	24
Einstellung während Zentrifugation.	33
Entsorgung.	48
Ersatzteile.	19
F	
Fehlermeldungen.	46
G	
Gerät	
desinfizieren.	44
reinigen.	43
Gummidichtung	
fetten.	45
I	
Integrale Zentrifugalbeschleunigung	
Integral RCF.	35
L	
Lid	
the lid.	24
Lieferumfang.	19
M	
Motorwelle	
fetten.	45
N	
NETZ-RESET.	47
Nicht vorgesehene Zweckbestimmung.	7
O	
Originalersatzteile.	19
P	
Personalunterweisung.	8
Personliche Schutzausrüstung.	7
Personnel qualifications.	7
Pflege	
Intervalle.	42
Programm	
ändern.	37
aufrufen.	36
eingeben.	37
laden.	36
Programmverknüpfung	
ändern.	41
erstellen.	40
löschen.	42
Zentrifugationslauf.	41
R	
Reinigung.	43
Reinigung und Desinfektion	
Hinweise.	43
Relative Zentrifugalbeschleunigung	
RCF.	35
Relative Zentrifugalbeschleunigung (RCF/RZB).	36
Rotor	
ausbauen.	25
beladen.	27, 28
einbauen.	25
Rotoren, Gehänge und Zubehör	
Berechnung der Anzahl der durchgeführten Laufzyklen.	46
Rotorerkennung.	37
Rücksendung.	19
Runtime.	34
S	
Schilder	
am Gerät.	14
auf der Verpackung.	14
Schleuderraum	
prüfen.	45
Schlüsselschalter.	33
Schutzausrüstung.	7
Sicherheitshinweise.	8
Storage conditions.	20
Störungsbehebung.	46
Symbole.	6
Systeminformationen	
abfragen.	38
T	
Tragzapfen	
fetten.	45
Transport condition.	20
Trouble shooting.	46
Typenschild.	13

V

Verantwortung des Betreibers.	8
Vorgesehene Zweckbestimmung.	6
Vorhersehbare Fehlanwendung.	7


W



Wartung.	45
Intervalle.	42

Z


Zentrifugation	
im Dauerlauf.	32
mit höherer Stoffdichte.	36
mit Zeitvorwahl.	33
Zentrifugationsdaten nach Einschalten.	39
Zentrifugiergefäße	
tauschen.	46
Zubehör.	19
desinfizieren.	44
mit begrenzter Verwendungsdauer.	46
prüfen.	45
reinigen.	44
Zwischenspeicher	
automatisch.	37

Rotoren und Zubehör / Rotors and accessories



4296		5051 + 5053							
Ausschwingrotor 6-fach / Swing out rotor 6-times									
  max. Laufzyklen / max. cycles 120000		 max. Beladung / max. load: 500 g							
		5262	5249	5243	5243 + 2x 6316	5242	5247	5227	5257
       		---	---	Falcon ®	---	---	---	---	
       		2)	2)	2)	3)	2)	2)	2)	2)
Kapazität / capacity	ml	100	100	50	50	25	7	5	6
Maße / dimensions	Ø x L mm	44 x 100	40 x 115	34 x 100	29 x 115	24 x 100	12 x 100	12 x 75	11 x 38
Anzahl p. Rotor / number p. rotor		6	6	12	12	30	120	120	240
Drehzahl / speed	RPM	4000							
RZB / RCF ²⁾		3291	3291	3291	3291	3291	3291	3309	2486 / 3363
Radius / radius	mm	184	184	184	184	184	184	185	139/188
<input checked="" type="checkbox"/> 9 (97%)	sec	33							
<input type="checkbox"/> 9	sec	50							
Temperatur / temperature	°C ¹⁾	0							

4296		5051 + 5053									
Ausschwingrotor 6-fach / Swing out rotor 6-times											
  max. Laufzyklen / max. cycles 120000		 max. Laufzyklen / max. cycles: 50000									
		5248-91	5247-91	5266	5258	5264	5227	5248			
       		10)	10)	3)							
       		2)	2)	2)							
Kapazität / capacity	ml	15	7	30	9 - 10	4 - 5,5	7,5 - 8,2	2,7 - 3	4,5 - 5		
Maße / dimensions	Ø / L mm ²	17 x 100	12 x 100	25 x 110	16 x 92	15 x 75	15 x 92	11 x 66	11 x 92		
Anzahl p. Rotor / number p. rotor		72	120	30	66	72		120	72		
Drehzahl / speed	RPM	4000									
RZB / RCF ²⁾		3291	3291	3291	3291	3309		3309	3291		
Radius / radius	mm	184	184	184	184	185		185	184		
<input checked="" type="checkbox"/> 9 (97%)	sec	33									
<input type="checkbox"/> 9	sec	50									
Temperatur / temperature	°C ¹⁾	0									

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
 2) Angaben des Röhrchenherstellers beachten.
 3) nicht mit Deckel 5053 verschließbar
 10) mit Dekantierhilfe
- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
 2) Observe the tube manufacturer's instructions.
 3) not possible to close the lid 5053
 10) with decanting aid

4296		5051 + 5053					
Ausschwingrotor 6-fach / Swing out rotor 6-times		 					
		max. Laufzyklen / max. cycles: 50000 max. Beladung / max. load: 500 g					



5259	6306	5248	5264	5267	5281		
 3)	 3)						
 ---		 ---					
 2)		 2)					
Kapazität / capacity ml		50	15	5 - 10	4 - 7	9	1,1 - 1,4
Maße / dimensions Ø x L mm		29 x 115	17 x 120	16 x 100	16 x 75	14 x 100	8 x 66
Anzahl p. Rotor / number p. rotor		12	42	72			120
Drehzahl / speed RPM					4000		
RZB / RCF ²⁾		3363	3434	3291	3309	3274	3363
Radius / radius mm		188	192	184	185	183	188
 9 (97%) sec					33		
 9 sec					50		
Temperatur / temperature °C ¹⁾					0		


4296		5051 + 5053					
Ausschwingrotor 6-fach / Swing out rotor 6-times		 					
		max. Laufzyklen / max. cycles: 50000 max. Beladung / max. load: 500 g					

5268							
							
 2)		 2)					
Kapazität / capacity ml		1 - 5	4 - 7	5	2,6 - 2,9	4,9	---
Maße / dimensions Ø x L mm		13 x 75	13 x 100	13 x 75	13 x 65	13 x 90	---
Anzahl p. Rotor / number p. rotor				72			---
Drehzahl / speed RPM				4000			---
RZB / RCF ²⁾				3345			---
Radius / radius mm				187			---
 9 (97%) sec				33			---
 9 sec				50			---
Temperatur / temperature °C ¹⁾				0			---












- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrchenherstellers beachten.
- 3) nicht mit Deckel 5053 verschließbar








- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 3) not possible to close the lid 5053

4296	5051	+	5280	+	5053		
Ausschwingrotor 6-fach / Swing out rotor 6-times  $\angle 90^\circ$ max. Laufzyklen / max. cycles 120000							
						+ max. Laufzyklen / max. cycles: 50000	max. Beladung / max. load: 500 g
						- - -	
	1662					1670	
							
	1663	1664	1665	1666	1667	1668	1663
							
Kapazität / capacity ml	1	2	4	8	3 x 2	4 x 1	1
Maße / dimensions Ø / A mm²	6,2 / 30	8,7 / 60	12,4 / 120	17,5 / 240	8,7 / 60	6,2 / 30	6,2 / 30
Anzahl p. Rotor / number p. rotor	12	12	12	12	12	12	12
Filterkarten / filter cards	1675	1675	1675	1676	1677	1678	1692
Drehzahl / speed RPM					4000		
RZB / RCF ²⁾					2290 / 3274		
Radius / radius mm					128 / 183		
<input checked="" type="checkbox"/> 9 (97%) sec					33		
<input type="checkbox"/> 9 sec					50		
Temperatur / temperature °C ¹⁾					0		

4296	5051	+	5280	+	5053		
Ausschwingrotor 6-fach / Swing out rotor 6-times  $\angle 90^\circ$ max. Laufzyklen / max. cycles 120000							
						+ max. Laufzyklen / max. cycles: 50000	max. Beladung / max. load: 500 g
						- - -	
	1670						
							
	1665	1666	1667	1668		---	---
						---	---
Kapazität / capacity ml	4	8	3 x 2	4 x 1		---	---
Maße / dimensions Ø / A mm²	12,4 / 120	17,5 / 240	8,7 / 60	6,2 / 30		---	---
Anzahl p. Rotor / number p. rotor	12	12	12	12		---	---
Filterkarten / filter cards	1692	1691	1694	1693		---	---
Drehzahl / speed RPM			4000			---	---
RZB / RCF ²⁾			2290 / 3274			---	---
Radius / radius mm			128 / 183			---	---
<input checked="" type="checkbox"/> 9 (97%) sec			33			---	---
<input type="checkbox"/> 9 sec			50			---	---
Temperatur / temperature °C ¹⁾			0			---	---




- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
 2) Angaben des Röhrchenherstellers beachten.
 4) Objekträger nur belastbar bis RZB 1100
- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
 2) Observe the tube manufacturer's instructions.
 4) Object slide will not stand RCF values exceeding 1100



4296		5092 + 5093							
Ausschwingrotor 6-fach / Swing out rotor 6-times									
  ↗ 90°		 +  max. Laufzyklen / max. cycles: 30000							
		max. Beladung / max. load: 500 g							
		mit Bioabdichtung / with bio-containment 5)							
5126	5125	5123	5124		5122		5128		
									
...
									
Kapazität / capacity ml	100	100	50	50	25	30	4	5	6
Maße / dimensions Ø x L mm	40 x 115	44 x 100	29 x 115	34 x 100	24 x 100	25 x 110	12 x 60	12 x 75	12 x 82
Anzahl p. Rotor / number p. rotor	6	6	12	6	24		72		72
Drehzahl / speed RPM					4000				
RZB / RCF ²⁾	3488	3488	3631	3488		3434			3542
Radius / radius mm	195	195	203	195		192			198
✓ 9 (97%) sec					33				
✗ 9 sec					50				
Temperatur / temperature °C ¹⁾					- 2				

4296		5092 + 5093							
Ausschwingrotor 6-fach / Swing out rotor 6-times									
  ↗ 90°		 +  max. Laufzyklen / max. cycles: 30000							
		max. Beladung / max. load: 500 g							
		mit Bioabdichtung / with bio-containment 5)							
5120		5121							
									
...	3)
									
Kapazität / capacity ml	7	4,5 - 5	4 - 7	2,6 - 2,9	9 - 10	10	5 - 10	8	15
Maße / dimensions Ø x L mm	12 x 100	11 x 92	13 x 100	13 x 65	16 x 92	15 x 102	16 x 100	16 x 125	17 x 100
Anzahl p. Rotor / number p. rotor		72				42			
Drehzahl / speed RPM					4000				
RZB / RCF ²⁾		3542				3542			
Radius / radius mm		198				198			
✓ 9 (97%) sec					33				
✗ 9 sec					50				
Temperatur / temperature °C ¹⁾					- 2				







- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrchenherstellers beachten.
- 3) nicht mit Deckel 5093 verschließbar
- 5) Nach DIN EN 61010, Teil 2 – 020. Die Hinweise für Bio-Sicherheitssysteme in den Kapiteln "Sicherheitshinweise" und "Pflege und Wartung" beachten.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 3) not possible to close the lid 5093
- 5) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".

4296	5092 + 5093
Ausschwingrotor 6-fach / Swing out rotor 6-times  ↙ 90° max. Laufzyklen / max. cycles 120000	 max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 500 g mit Bioabdichtung / with bio-containment 5) 5136  
Kapazität / capacity ml	10 15 --- 4 - 5,5 7,5 - 8,2 9 - 10 10 4 - 7 5 - 10
Maße / dimensions Ø x L mm	16 x 80 17 x 100 --- 15 x 75 15 x 92 16 x 92 15 x 102 16 x 75 16 x 100
Anzahl p. Rotor / number p. rotor	48
Drehzahl / speed RPM	4000
RZB / RCF ²⁾	3488
Radius / radius mm	195
↙ 9 (97%) sec	33
↖ 9 sec	50
Temperatur / temperature °C ¹⁾	- 2





4296	5092 + 5093
Ausschwingrotor 6-fach / Swing out rotor 6-times  ↙ 90° max. Laufzyklen / max. cycles 120000	 max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 500 g mit Bioabdichtung / with bio-containment 5) 5137  
Kapazität / capacity ml	1 - 5 4 - 7 4,9 1,1 - 1,4 2,6 - 2,9 2,7 - 3 4,5 - 5 5
Maße / dimensions Ø x L mm	13 x 75 13 x 100 13 x 90 8 x 66 13 x 65 11 x 66 11 x 92 13 x 75
Anzahl p. Rotor / number p. rotor	48
Drehzahl / speed RPM	4000
RZB / RCF ²⁾	3488
Radius / radius mm	195
↙ 9 (97%) sec	33
↖ 9 sec	50
Temperatur / temperature °C ¹⁾	- 2









- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrchenherstellers beachten.
- 5) Nach DIN EN 61010, Teil 2 – 020. Die Hinweise für Bio-Sicherheitssysteme in den Kapiteln "Sicherheitshinweise" und "Pflege und Wartung" beachten.
- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 5) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".

4296	5092 + 5093								5092								
Ausschwingrotor 6-fach / Swing out rotor 6-times  max. Laufzyklen / max. cycles 120000	 max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 500 g mit Bioabdichtung / with bio-containment 5)																
1791 5134 5135 5129 5138 6319 6319																	
 ---																	
 2) ---																	
Kapazität / capacity ml	250	25	50	15	1.1 – 1.4	2,7 - 3	2,6 – 2,9	1 – 5	250	290							
Maße / dimensions Ø x L mm	65 x 115	25 x 90	29 x 115	17 x 120	8 x 66	11 x 66	13 x 65	13 x 75	62 x 122	62 x 137							
Anzahl p. Rotor / number p. rotor	6	18	12	42			72		6	6							
Drehzahl / speed RPM					4000												
RZB / RCF ²⁾	3631	3363	3560	3631			3077		3631	3631							
Radius / radius mm	203	188	199	203			172		203	203							
<input checked="" type="checkbox"/> 9 (97%) sec							33										
<input checked="" type="checkbox"/> 9 sec							50										
Temperatur / temperature °C ¹⁾							- 2										
4294	4290 + 4291																
Ausschwingrotor 4-fach / Swing out rotor 4-times  max. Laufzyklen / max. cycles 40000	 max. Laufzyklen / max. cycles: 30000 (4500 - 4001 RPM) 45000 (4000 – 3501 RPM) 60000 (3500 - 50 RPM)																
4273																	
 ---																	
 2) ---																	
Kapazität / capacity ml	5	6	7	2,6 - 2,9	4,9	1 - 5	4 - 7		---								
Maße / dimensions Ø x L mm	12 x 75	12 x 82	12 x 100	13 x 65	13 x 90	13 x 75	13 x 100		---								
Anzahl p. Rotor / number p. rotor				200					---								
Drehzahl / speed RPM					4500												
RZB / RCF ²⁾					4551												
Radius / radius mm					201												
<input checked="" type="checkbox"/> 9 (97%) sec					115												
<input checked="" type="checkbox"/> 9 sec					116												
Temperatur / temperature °C ¹⁾					2												





- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrchenherstellers beachten.
- 5) Nach DIN EN 61010, Teil 2 – 020. Die Hinweise für Bio-Sicherheitssysteme in den Kapiteln "Sicherheitshinweise" und "Pflege und Wartung" beachten.
- 12) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.








- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 5) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".
- 12) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.

4294		4290 + 4291						
Ausschwingrotor 4-fach / Swing out rotor 4-times		 max. Laufzyklen / max. cycles: 30000 (4500 - 4001 RPM) 45000 (4000 - 3501 RPM) 60000 (3500 - 50 RPM)						
 ↗ 90°		max. Beladung / max. load: 1200 g mit Bioabdichtung / with bio-containment 5)						
max. Laufzyklen / max. cycles 40000		4310	4311	4313	4314	4321	---	---
							---	---
							---	---
Kapazität / capacity ml		12	10	9 - 10	50	15	50	---
Maße / dimensions Ø x L mm		16,8 x 100	15 x 102	16 x 92	29 x 115	17 x 120	29 x 115	---
Anzahl p. Rotor / number p. rotor		112	132		32	68	32	---
Drehzahl / speed RPM		4500						
RZB / RCF ²⁾		4437	4573	4528	4618	4618	---	---
Radius / radius mm		196	202	200	204	204	---	---
 9 (97%) sec		115						
 9 sec		116						
Temperatur / temperature °C ¹⁾		2						

4294		4290 + 4291						
Ausschwingrotor 4-fach / Swing out rotor 4-times		 max. Laufzyklen / max. cycles: 30000 (4500 - 4001 RPM) 45000 (4000 - 3501 RPM) 60000 (3500 - 50 RPM)						
 ↗ 90°		max. Beladung / max. load: 1200 g mit Bioabdichtung / with bio-containment 5)						
max. Laufzyklen / max. cycles 40000								
								
Kapazität / capacity ml		9	14	15	4 - 5,5	7,5 - 8,2	4 - 7	5 - 10
Maße / dimensions Ø x L mm		14 x 100	16 x 101	17 x 100	15 x 75	15 x 92	16 x 75	16 x 100
Anzahl p. Rotor / number p. rotor		168						
Drehzahl / speed RPM		4500						
RZB / RCF ²⁾		4551						
Radius / radius mm		201						
 9 (97%) sec		115						
 9 sec		116						
Temperatur / temperature °C ¹⁾		2						


- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrchenherstellers beachten.
- 5) Nach DIN EN 61010, Teil 2 – 020. Die Hinweise für Bio-Sicherheitssysteme in den Kapiteln "Sicherheitshinweise" und "Pflege und Wartung" beachten.
- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 5) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".


4294	4290							
Ausschwingrotor 4-fach / Swing out rotor 4-times								
	max. Laufzyklen / max. cycles: 30000 (4500 - 4001 RPM) 45000 (4000 - 3501 RPM) 60000 (3500 - 50 RPM)							
max. Laufzyklen / max. cycles 40000	max. Beladung / max. load: 1200 g							
4339	4323	4320	---	---	---	---	---	---
			---	---	---	---	---	---
---	---	---	---	---	---	---	---	---
			---	---	---	---	---	---
Kapazität / capacity ml	50	50	15	---	---	---	---	---
Maße / dimensions Ø x L mm	29 x 115	29 x 115	17 x 120	---	---	---	---	---
Anzahl p. Rotor / number p. rotor	40	40	112	---	---	---	---	---
Drehzahl / speed RPM	4500							
RZB / RCF ²⁾	4528	4618	4618	---	---	---	---	---
Radius / radius mm	200	204	204	---	---	---	---	---
<input checked="" type="checkbox"/> 9 (97%) sec	115							
<input checked="" type="checkbox"/> 9 sec	116							
Temperatur / temperature °C ¹⁾	2							

4294	4295-A + 4229-B								
Ausschwingrotor 4-fach / Swing out rotor 4-times									
	max. Laufzyklen / max. cycles: 15000 (4500 - 4001 RPM) 50000 (4000 - 50 RPM)								
max. Laufzyklen / max. cycles 40000	max. Beladung / max. load: 1060 g								
4226	4225	4224	4241	4245-A	4213				
					3)				
---	---	---	---	---	---	---	---	---	
					2)		2)	2)	
Kapazität / capacity ml	0.8	1.5	2.0	4	25	50	6	7	4.5 - 5
Maße / dimensions Ø x L mm	8 x 45	11 x 38	10 x 88	25 x 90	29 x 115	12 x 82	12 x 100	11 x 92	
Anzahl p. Rotor / number p. rotor	312	336	252	32	32			192	
Drehzahl / speed RPM	4500								
RZB / RCF ²⁾	4777	3690 / 4887	4777	4777	4958			4777	
Radius / radius mm	211	163 / 215	211	211	219			211	
<input checked="" type="checkbox"/> 9 (97%) sec	115								
<input checked="" type="checkbox"/> 9 sec	116								
Temperatur / temperature °C ¹⁾	6								

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
 2) Angaben des Röhrchenherstellers beachten.
 3) nicht mit Deckel 4229-B verschließbar



- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
 2) Observe the tube manufacturer's instructions.
 3) not possible to close the lid 4229-B


4294	4295-A + 4229-B												
Ausschwingrotor 4-fach / Swing out rotor 4-times													
													
max. Laufzyklen / max. cycles 40000													
90°													
Kapazität / capacity ml	5	6	2,7 - 3	15	7,5 - 8,2	5 - 10	10	4 - 7					
Maße / dimensions Ø x L mm	12 x 75	12 x 82	11 x 66	17 x 100	15 x 92	16 x 100	15 x 102	16 x 75					
Anzahl p. Rotor / number p. rotor	192			120			120						
Drehzahl / speed RPM	4500												
RZB / RCF ²⁾	4777			4777			4777						
Radius / radius mm	211			211			211						
9 (97%) sec	115												
9 sec	116												
Temperatur / temperature °C ¹⁾	6												

4294	4295-A + 4229-B							
Ausschwingrotor 4-fach / Swing out rotor 4-times								
								
max. Laufzyklen / max. cycles 40000								
90°								
Kapazität / capacity ml	4 - 5,5	---	50	100	94	250	250	---
Maße / dimensions Ø x L mm	15 x 75	---	34 x 100	40 x 115	38 x 102	62 x 122	56 x 144	---
Anzahl p. Rotor / number p. rotor	120	---	24	16	4	4	4	---
Drehzahl / speed RPM	4500							
RZB / RCF ²⁾	4777	---	4777	4777	4777	4641	4641	---
Radius / radius mm	211	---	211	211	211	205	205	---
9 (97%) sec	115							
9 sec	116							
Temperatur / temperature °C ¹⁾	6							

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrchenherstellers beachten.
- 12) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 12) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.

4294	4295-A + 4229-B							
Ausschwingrotor 4-fach / Swing out rotor 4-times								
								
max. Laufzyklen / max. cycles: 40000								
90°	4220	4222	4223	---	---	---	---	---
			---	---	---	---	---	---
---	---	---	---	---	---	---	---	---
			2)	2)	2)	2)	2)	---
Kapazität / capacity ml	9 - 10	12	4 - 7	7	9	12	8	---
Maße / dimensions Ø x L mm	16 x 92	16,8 x 100	13 x 100	12 x 100	14 x 100	16 x 101	16 x 125	---
Anzahl p. Rotor / number p. rotor	64	120	120	120	100	100	100	---
Drehzahl / speed RPM	4500							---
RZB / RCF ²⁾	4777	4777	4777	4777	4777	4777	4777	---
Radius / radius mm	211	211	211	211	211	211	211	---
<input checked="" type="checkbox"/> 9 (97%) sec	115							---
<input type="checkbox"/> 9 sec	116							---
Temperatur / temperature °C ¹⁾	6							---

4294	4295-A + 4229-B							
Ausschwingrotor 4-fach / Swing out rotor 4-times								
								
90°								
max. Laufzyklen / max. cycles: 40000								
	4249	4222-93	4258	---	---	---	---	---
			---	---	---	---	---	---
---	---	---	---	---	---	---	---	---
			12)	12)	12)	12)	12)	12)
Kapazität / capacity ml	---	50	2,6 – 2,9	1 - 5	750	750	500	650
Maße / dimensions Ø x L mm	---	29 x 115	13 x 65	13 x 75	96 x 135	97 x 152	96 x 147	97 x 139
Anzahl p. Rotor / number p. rotor	---	24	120	120	120	120	4	4
Drehzahl / speed RPM	4500							---
RZB / RCF ²⁾	---	4867	4777	4777	4777	4777	4958	4958
Radius / radius mm	---	215	211	211	211	211	219	219
<input checked="" type="checkbox"/> 9 (97%) sec	115							---
<input type="checkbox"/> 9 sec	116							---
Temperatur / temperature °C ¹⁾	6							---

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrchenherstellers beachten.
- 12) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 12) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.

4294	4295-A + 4229-B				4295-A
Ausschwingrotor 4-fach / Swing out rotor 4-times					
	max. Laufzyklen / max. cycles: 15000 (4500 - 4001 RPM) 50000 (4000 - 50 RPM)				
↙ 90°	max. Beladung / max. load: 1060 g				
max. Laufzyklen / max. cycles 40000	6322	---	4232	4215	SK15.16





	(3)		(3)	--- / 3)	(3)
Corning	---	---	---	---	---

					Blut-Dopingtest
Kapazität / capacity ml	250	---	15	25	---
Maße / dimensions Ø x L mm	60 x 162	---	17 x 120	24 x 100	25 x 110
Anzahl p. Rotor / number p. rotor	4	---	92	44	24
Drehzahl / speed RPM	4500				
RZB / RCF ²⁾	4777	---	4958	4777	---
Radius / radius mm	211	---	219	211	---
↙ 9 (97%) sec	115				
↙ 9 sec	116				
Temperatur / temperature °C ¹⁾	6				

4294	4298-A				4293				
Ausschwingrotor 4-fach / Swing out rotor 4-times									
	max. Laufzyklen / max. cycles 50000 max. Beladung / max. load: 1150 g								
↙ 90°	max. Laufzyklen / max. cycles 50000 max. Beladung / max. load: 1150 g								
max. Laufzyklen / max. cycles 40000	---	4237-A	4244-A						

1-fach/ 1-times	4-fach/ 4-times	4-fach/ 4-times	1-fach/ 1-times	3-fach/ 3-times	2-fach/ 2-times				
Kapazität / capacity ml	1000	450	500	750	450				
Maße / dimensions Ø x L mm	---	---	---	---	---				
Anzahl p. Rotor / number p. rotor	4	4	4	4	4				
Drehzahl / speed RPM	4500								
RZB / RCF ²⁾	5071	5003		5252					
Radius / radius mm	224	221		232					
↙ 9 (97%) sec	115								
↙ 9 sec	116								
Temperatur / temperature °C ¹⁾	3			2					

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
 2) Angaben des Röhrchenherstellers beachten.
 3) nicht mit Deckel 4229-B verschließbar
- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
 2) Observe the tube manufacturer's instructions.
 3) not possible to close the lid 4229-B

4294		4280 + 5629								
Ausschwingrotor 4-fach / Swing out rotor 4-times		 max. Laufzyklen / max. cycles: 50000 max. Beladung / max. load: 690 g mit Bioabdichtung / with bio-containment 5)								
 max. Laufzyklen / max. cycles: 40000										
MTP 		4279 								
MTP 		CP	MS	DWP	QP	96-PCR-Platte / plate	PCR-Strips			
										
Kapazität / capacity ml		---	---	---	---	---	---	0,2		
Maße / dimensions mm		86x128x15	86x128x17,5	86x128x22	86x128x46	86x128x44,5	86x128x 83	82x124x20		
Anzahl p. Rotor / number p. rotor		24	24	20	4	8	4	48 x 8		
Drehzahl / speed RPM		4500								
RZB / RCF ²⁾		4573								
Radius / radius mm		202								
 9 (97%) sec		115								
 9 sec		116								
Temperatur / temperature °C ¹⁾		5								
4294		4257	4254 + 4255 / 4255-P ⁸⁾							
Ausschwingrotor 4-fach / Swing out rotor 4-times			 + 9)							
 max. Laufzyklen / max. cycles: 40000		max. Laufzyklen / max. cycles: 20000 max.							max. Laufzyklen / max. cycles: 30000	
 max. Laufzyklen / max. cycles: 40000		max. Beladung / max. load: 800 g (4500 – 4021 RPM) 1000 g (4020 – 3671 RPM) 1200 g (3670 - 50 RPM)								
...		...								
...		4259-A	---	---	---	4449	4430			
...			---	---	---					
Hitachi-Racks		0554	0512	4239	Corning	Corning	Nagene	Nunc		
		 12)	 12)	 12)						
Kapazität / capacity ml		650	750	1000	500	250	175	200		
Maße / dimensions Ø x L mm		20 x 118 x 70	20 x 118 x 70	97 x 139	97 x 152	96 x 176	96 x 147	60 x 162	61,5 x 139,2	
Anzahl p. Rotor / number p. rotor		20	20	4	4	4	4	4	4	
Drehzahl / speed RPM		4500								
RZB / RCF ²⁾		4822	4867	5184						
Radius / radius mm		213	215	229						
 9 (97%) sec		115								
 9 sec		116								
Temperatur / temperature °C ¹⁾		7	6							

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrchenherstellers beachten.
- 5) Nach DIN EN 61010, Teil 2 – 020. Die Hinweise für Bio-Sicherheitssysteme in den Kapiteln "Sicherheitshinweise" und "Pflege und Wartung" beachten.
- 7) Bei Verwendung von Entnahmehilfe 4259-A die Einlage (E2435) aus Gehänge 4257 herausnehmen. Entnahmehilfe 4259-A nur voll beladen zentrifugieren.
- 8) 4255-P: spezielle Oberflächenbehandlung für höchste hygienische Ansprüche
- 9) ohne Deckel
- 12) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.
- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 5) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".
- 7) When using removal frame 4259-A please take insert (E2435) out of hanger 4257. Centrifuge removal frame 4259-A only when fully loaded.
- 8) 4255-P: special surface treatment for highest hygienic requirements
- 9) without lid
- 12) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.



MTP Mikrotiterplatte / Microtitre plate


CP Kulturplatte / Culture plate

DWP Deep Well Platte / Deep well plate

MS Micronic System / Micronic system


QP Filterplatte / Filter plate

4294	4254	4254 + 4255 / 4255-P ⁸⁾	
 <p>Ausschwingrotor 4-fach / Swing out rotor 4-times max. Laufzyklen / max. cycles 40000</p>		 <p>max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 800 g (4500 – 4021 RPM) 1000 g (4020 – 3671 RPM) 1200 g (3670 - 50 RPM)</p>	

---		4432	4433
---			
4255 / 4255-P ⁸⁾		4434	
---			
Kapazität / capacity ml	1000	1,5	2,0
Maße / dimensions Ø x L mm	98 x 138	11 x 38	12 x 75
Anzahl p. Rotor / number p. rotor	4	168	120
Drehzahl / speed RPM			4500
RZB / RCF ²⁾	5184	3600/4686	4618
Radius / radius mm	229	159/207	204
9 (97%) sec			115
9 sec			116
Temperatur / temperature °C ¹⁾			6
4294	4254	4254 + 4255 / 4255-P ⁸⁾	
 <p>Ausschwingrotor 4-fach / Swing out rotor 4-times max. Laufzyklen / max. cycles 40000</p>		 <p>max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 800 g (4500 – 4021 RPM) 1000 g (4020 – 3671 RPM) 1200 g (3670 - 50 RPM)</p>	


---		4434	
---			
Kapazität / capacity ml	15	10	8
Maße / dimensions Ø x L mm	17 x 100	16 x 80	16 x 81
Anzahl p. Rotor / number p. rotor			76
Drehzahl / speed RPM			4500
RZB / RCF ²⁾			4777
Radius / radius mm			211
9 (97%) sec			115
9 sec			116
Temperatur / temperature °C ¹⁾			6

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrchenherstellers beachten.
- 3) 4255 nicht mit Deckel verschließbar
- 8) 4255-P: spezielle Oberflächenbehandlung für höchste hygienische Ansprüche
- 13) Maximale Beladung 800g. Bei einer Beladung über 800g muss die Drehzahl reduziert werden, siehe Beschriftung auf dem Becher. Berechnung der reduzierten Drehzahl siehe Kapitel "Zentrifugation von Stoffen oder Stoffgemischen mit einer höheren Dichte als 1,2 kg/dm3".
- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 3) 4255 not possible to close the lid
- 8) 4255-P: special surface treatment for highest hygienic requirements
- 13) Maximum load 800g. With a load higher than 800g the speed has to be reduced, see label on the bucket. Calculation of the reduced speed see chapter "Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm3".




















4294	4254 + 4255 / 4255-P ⁸⁾													
Ausschwingrotor 4-fach / Swing out rotor 4-times	  max. Laufzyklen / max. cycles: 30000													
	max. Beladung / max. load: 800 g (4500 – 4021 RPM) 1000 g (4020 – 3671 RPM) 1200 g (3670 - 50 RPM)													

	4435		4437	4438		4438 + 0726								
														

Kapazität / capacity ml	2,6 – 2,9	4,9	1 – 5	4 – 7	15	25	30	25						
Maße / dimensions Ø x L mm	13 x 65	13 x 90	13 x 75	13 x 100	17 x 120	25 x 90	25 x 110	24 x 100						
Anzahl p. Rotor / number p. rotor	84		48		28		28							
Drehzahl / speed RPM	4500													
RZB / RCF ²⁾	4618		4890		4709		4505							
Radius / radius mm	204													
<input checked="" type="checkbox"/> 9 (97%) sec	115													
<input type="checkbox"/> 9 sec	116													
Temperatur / temperature °C ¹⁾	6													
4294	4254 + 4255 / 4255-P ⁸⁾													
Ausschwingrotor 4-fach / Swing out rotor 4-times	  max. Laufzyklen / max. cycles: 30000													
	max. Beladung / max. load: 800 g (4500 – 4021 RPM) 1000 g (4020 – 3671 RPM) 1200 g (3670 - 50 RPM)													
	4439		4440	4441	4442	4443		---						
								---						
								---						
														
Kapazität / capacity ml	50	225	175	50	100	250	290	---						
Maße / dimensions Ø x L mm	34 x 100	61 x 137	61 x 118	29 x 115	44 x 100	62 x 122	62 x 137	---						
Anzahl p. Rotor / number p. rotor	16	4	20	8	4		---							
Drehzahl / speed RPM	4500													
RZB / RCF ²⁾	4573	5184	4890		4551	5003		---						
Radius / radius mm	202	229	216		201	221		---						
<input checked="" type="checkbox"/> 9 (97%) sec	115													
<input type="checkbox"/> 9 sec	116													
Temperatur / temperature °C ¹⁾	6													

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrchenherstellers beachten.
- 3) 4255 nicht mit Deckel verschließbar
- 8) 4255-P: spezielle Oberflächenbehandlung für höchste hygienische Ansprüche
- 12) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.
- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 3) 4255 not possible to close the lid
- 8) 4255-P: special surface treatment for highest hygienic requirements
- 12) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.

4282		4285-A								
Ausschwingrotor 2-fach / Swing out rotor 2-times   		 max. Laufzyklen / max. cycles: 12000 (3600 - 3001 RPM) 30000 (3000 - 50 RPM) max. Beladung / max. load: 2320 g ---  4281 + 2x 1485								
  		 max. Laufzyklen / max. cycles: 100000 --- MTP MTP MS CP DWP Microtest-platten / plate Terasaki 96-PCR-Platte / plate PCR-Strips								
  		 max. Laufzyklen / max. cycles: 100000 --- Kapazität / capacity ml --- --- --- --- --- --- --- 0,2 Maße / dimensions mm 86x128x15 86x128x17,5 86x128x46 86x128x22 86x128x44,5 59x84x11 82x124x20 --- Anzahl p. Rotor / number p. rotor 16 16 4 12 4 4 4 48 x 4 Drehzahl / speed RPM 3600 RZB / RCF 2) 2434 Radius / radius mm 168 <input checked="" type="checkbox"/> 9 (97%) sec 87 <input checked="" type="checkbox"/> 9 sec 94 Temperatur / temperature °C 1) - 5								
  		 max. Laufzyklen / max. cycles: 12000 (3600 - 3001 RPM) 30000 (3000 - 50 RPM) max. Beladung / max. load: 2320 g ---  4285-A								
  		 max. Laufzyklen / max. cycles: 12000 (3600 - 3001 RPM) 30000 (3000 - 50 RPM) max. Beladung / max. load: 2320 g --- 4263-A SK 01.14 4283-B 4287-B 4288-A SK 25.10 + SK 25.10-1 SK 06.21-01 + SK 06.21-02 SK 32.07        S-Monovette® / Rack Sarstedt AutoMate™ Rack / Beckman Coulter Olympus-Racks Hitachi-Racks Behring Rack --- --- --- Positionen/positions 50 50 --- --- --- --- --- --- Kapazität / capacity ml --- --- --- --- --- --- --- --- Maße / dimensions mm 209x109x45 209x109x45 20x41x176 20x70x118 25x60x193 --- --- --- Anzahl p. Rotor / number p. rotor 2 2 12 20 10 2 2 10 Drehzahl / speed RPM 3600 RZB / RCF 2) 2579 2579 2652 2652 2652 2594 2492 2652 Radius / radius mm 178 178 183 183 183 179 172 183 <input checked="" type="checkbox"/> 9 (97%) sec 87 <input checked="" type="checkbox"/> 9 sec 94 Temperatur / temperature °C 1) - 5								

1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)

2) Angaben des Röhrchenherstellers beachten.

1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)

2) Observe the tube manufacturer's instructions.


MTP Mikrotiterplatte /
Microtitre plate


CP Kulturplatte /
Culture plate

DWP Deep Well Platte /
Deep well plate



MS Micronic System /
Micronic system

QP Filterplatte /
Filter plate

<p>4246</p> <p>Winkelrotor 6-fach / Angle rotor 6-times</p>  <p>max. Laufzyklen / max. cycles 15000</p> <p>Hülsen / reduction (6x inclusive) max. Laufzyklen: 15000 einsetzbar bis / usable until: 5 Jahre / year mit Bioabdichtung / with bio-containment 5)</p> <p>Kapazität / capacity ml 94 85 50 25 30 7,5 - 8,2 9 - 10 10</p> <p>Maße / dimensions Ø x L mm 38 x 110 38 x 106 29 x 107 24 x 100 26 x 95 15 x 92 16 x 92 15 x 102</p> <p>Anzahl p. Rotor / number p. rotor 6 6 6 6 6 6 6 6</p> <p>Drehzahl / speed RPM 11500</p> <p>RZB / RCF 2) 18038 17299 16560 17003</p> <p>Radius / radius mm 122 117 112 115</p> <p>✓ 9 (97%) sec 64</p> <p>✗ 9 sec 64</p> <p>Temperatur / temperature °C 16) 2</p>	<p>---</p> <p>---</p> <p>---</p> <p>---</p> <p>---</p> <p>---</p> <p>---</p> <p>1446 1447 1451</p> <p>---</p>
---	---

<p>4246</p> <p>Winkelrotor 6-fach / Angle rotor 6-times</p>  <p>max. Laufzyklen / max. cycles 15000</p> <p>Hülsen / reduction (6x inclusive) max. Laufzyklen: 15000 einsetzbar bis / usable until: 5 Jahre / year mit Bioabdichtung / with bio-containment 5)</p> <p>Kapazität / capacity ml 5 - 10 15 15 50 75 94 10 15</p> <p>Maße / dimensions Ø x L mm 16 x 100 17 x 100 17 x 100 34 x 100 35 x 105 38 x 102 16 x 80 17 x 120</p> <p>Anzahl p. Rotor / number p. rotor 6 6 6 6 6 6 12 6</p> <p>Drehzahl / speed RPM 11500</p> <p>RZB / RCF 2) 17003 17743 18038 17003 17299</p> <p>Radius / radius mm 115 120 122 115 117</p> <p>✓ 9 (97%) sec 64</p> <p>✗ 9 sec 64</p> <p>Temperatur / temperature °C 16) 2</p>	<p>---</p> <p>---</p> <p>---</p> <p>---</p> <p>---</p> <p>1451 1463 1448 1466</p> <p>---</p>
---	--

- 2) Angaben des Röhrchenherstellers beachten.
- 5) Nach DIN EN 61010, Teil 2 – 020. Die Hinweise für Bio-Sicherheitssysteme in den Kapiteln "Sicherheitshinweise" und "Pflege und Wartung" beachten.
- 14) Maße mit Deckel 38 x 110 mm
- 16) Niedrigste Probentemperatur bei Vorkühlung und maximaler Drehzahl
- 2) Observe the tube manufacturer's instructions.
- 5) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".
- 14) Dimensions with cap 38 x 110 mm
- 16) Lowest sample temperature with pre-cooling and maximum speed (only with cooling centrifuges)

4246	---							
Winkelrotor 6-fach / Angle rotor 6-times  max. Laufzyklen / max. cycles 15000 Hülsen / reduction (6x inclusive) max. Laufzyklen: 15000 einsetzbar bis / usable until: 5 Jahre / year mit Bioabdichtung / with bio-containment ⁵⁾ Kapazität / capacity ml 1,5 2,0 0,5 3 4 5 50 --- Maße / dimensions Ø x L mm 11 x 38 11 x 38 10,7 x 44,5 10 x 60 12 x 40 17 x 59 29 x 115 --- Anzahl p. Rotor / number p. rotor 24 24 6 6 --- Drehzahl / speed RPM 11500 RZB / RCF ²⁾ 17299 17003 16856 17595 --- Radius / radius mm 117 115 114 119 --- <input checked="" type="checkbox"/> 9 (97%) sec <input type="checkbox"/> 9 sec Temperatur / temperature °C ¹⁶⁾ 2 --- --- --- 	--- --- --- 1449 1403 1476 1454 ---                							
4266 Winkelrotor 6-fach / Angle rotor 6-times  max. Laufzyklen / max. cycles 15000 Hülsen / reduction (6x inclusive) max. Laufzyklen: 15000 einsetzbar bis / usable until: 5 Jahre / year mit Bioabdichtung / with bio-containment ⁵⁾ Kapazität / capacity ml 250 10 30 25 50 94 85 94 15 Maße / dimensions Ø x L mm 61 x 122 16 x 80 26 x 95 24 x 100 29 x 107 38 x 110 38 x 106 38 x 102 17 x 100 Anzahl p. Rotor / number p. rotor 6 48 18 6 6 6 6 42 Drehzahl / speed RPM 9500 RZB / RCF ²⁾ 14025 13420 12915 12108 12310 13319 Radius / radius mm 139 133 128 120 122 132 <input checked="" type="checkbox"/> 9 (97%) sec <input type="checkbox"/> 9 sec Temperatur / temperature °C ¹⁶⁾ 2 --- --- --- 	--- --- --- 5641 5642 5643 5644 5646           <img alt="Diagram of a 6-place angle rotor with a 25° angle." data-b							