

SMALL CENTRIFUGE

MANUAL CENTRIFUGE

OPERATE INDEPENDENTLY ANYWHERE

The adaptable manual centrifuge can be simply clamped to any benchtop, with no need for electricity. Its self-lubricating drive system makes it virtually maintenance-free.



HIGHLIGHTS

- max. RPM: 3,000 min⁻¹
- max. RCF: 1,298
- max. capacity: 4 x 15 ml
- self-contained
- maintenance-free due to its self-lubricating drive system
- clamping device for simple table mounting
- choice of 2 rotors

TECHNICAL DATA

	MANUAL CENTRIFUGE without rotor	Swing-out rotor, 4-place	Swing-out rotor, 4-place
angle	-	90°	50°
max. capacity	-	4 x 15 ml (conical only)	4 x 15 ml
max. RPM	-	3,000 min ⁻¹	3,000 min ⁻¹
max. RCF	-	1,298	1,077
dimensions (WxDxH)	140 x 175 x 285 mm	-	-
weight	approx. 0.9 kg	-	-
Cat. No.	1011	1014	1025

COMPANY CERTIFICATIONS

