

SMALL CENTRIFUGE

# MANUAL CENTRIFUGE

## OPERATE INDEPENDENTLY ANYWHERE

The adaptable manual centrifuge can be simply clamped to any benchtop, with no need for electricity. Its self-lubricating drive system makes it virtually maintenance-free.



## HIGHLIGHTS

- max. RPM: 3,000 min<sup>-1</sup>
- max. RCF: 1,298
- max. capacity: 4 x 15 ml
- self-contained
- maintenance-free due to its self-lubricating drive system
- clamping device for simple table mounting
- choice of 2 rotors

## TECHNICAL DATA

	MANUAL CENTRIFUGE without rotor	Rotor 1014 4-place	Rotor 1025 4-place
angle	-	90°	50°
max. capacity	-	4 x 15 ml (conical only)	4 x 15 ml
max. RPM	-	3,000 min <sup>-1</sup>	3,000 min <sup>-1</sup>
max. RCF	-	1,298	1,077
dimensions (WxDxH)	140 x 175 x 285 mm	-	-
weight	approx. 0.9 kg	-	-
<b>Cat. No.</b>	<b>1011</b>	<b>1014</b>	<b>1025</b>

## CERTIFICATES / REGISTRATIONS

