

Benchtop centrifuges

ROTOFIX 46 | 46 H

Robust industrial unit

These are ideal centrifuges for research and industrial laboratories. They are robust and can withstand intensive usage, even in the field. Specialty buckets are able to hold ASTM pear-shaped and conical tubes for petroleum testing.

Heated option provides chamber temperature up to +90 °C (ROTOFIX 46 H).



■ Highlights

- max. RPM: 4,000 min⁻¹ | 2.000 min⁻¹
Adjustable in increments of 10
- max. RCF: 3,095 | 984
- max. capacity: 4 x 290 ml | 4 x 100 ml
- robust and versatile benchtop centrifuge
- choice of 4 rotors
- maximum noise level of ≤ 58 dB(A)
- impulse button for short centrifugation
- easy operation with keypad
- Model 46 H heatable from +10 to +90 °C

■ Features

- metal housing and lid
- viewing port in the lid
- twist lock
- lid dropping protection
- emergency lid lock release
- stainless steel chamber
- automatic rotor recognition
- brushless drive
- error display
- imbalance switch-off
- control panel with actual values of all parameters
- auditory message after completion of the centrifugation run

Swing-out rotor, 4-place | 4474


Rotor

max. RPM max. RCF	2,000 min ⁻¹ 984
max. capacity	4 x 100 ml
run-up run-down, braked in sec	16 40
angle max. noise level	90° 58 dB (A)
Cat. No.	4474




Bucket

Cat. No.	4275
----------	------



Tubes

	ASTM ^{2,1)}	0528 ^{2,1)}	ASTM ^{2,1)}	0531 ^{2,1)}	Babcock
capacity in ml	50	100	100	100	50
Ø x L in mm	45 x 130	58 x 161	44 x 168	37 x 200	36.5 x 185
max. RCF ^{2,1)}	984	984	961	961	912
radius in mm	220	220	215	215	204



Adapter

boring Ø x L in mm	-	-	-	-	-
tubes per rotor	4	4	4	4	4
Cat. No.	4278-A	0771	4277	4276-B	0703-A



Swing-out rotor, 4-place | 5694 (not for 46 H)



Rotor

max. RPM max. RCF	4,000 min ⁻¹ 3,095
max. capacity	4 x 100 ml
run-up run-down, braked in sec	40 45
angle max. noise level	90° 58 dB (A)
Cat. No.	5694




Bucket

lid	5053
Cat. No.	5051




Tubes

	tubes ²⁾										tubes with screw cap			
capacity in ml	3	5	6	7	9	15	25	50	100	100	15	50	12	50
Ø x L in mm	10 x 60	12 x 75	12 x 82	12 x 100	14 x 100	17 x 100	24 x 100	34 x 100	44 x 100	40 x 115	17 x 120	29 x 115	17 x 100	29 x 115
max. RCF ²⁾	2,737	2,773	2,773	2,755	2,773	2,755	2,755	2,755	2,755	2,755	2,898	2,844	2,898	2,755
radius in mm	153	155	155	154	155	154	154	154	154	154	162	158	162	154



Adapter

boring Ø x L in mm	11 x 44	12.5 x 39.5	12.5 x 39.5	12.5 x 82.5	16 x 50	17.5 x 82.5	26 x 82.5	36 x 82.5	45.5 x 82.5	42 x 82.5	17 x 90	30 x 82.5	17 x 90	36 x 82.5
tubes per rotor	80	80	80	80	48	48	20	8	4	4	28	8	28	8
Cat. No.	5267	5227	5227	5247¹⁵⁾	5264	5248¹⁶⁾	5242	5243	5262	5249	6306	5259	6306	5243⁴⁾



Swing-out rotor, 4-place | 5694 (not for 46 H)


Rotor

max. RPM max. RCF ²⁾	4,000 min ⁻¹ 3,095
max. capacity	4 x 290 ml
run-up run-down, braked in sec	40 45
angle max. noise level	90° 58 dB (A)
Cat. No.	5694



Bucket

lid bioseal ³⁾	5093
Cat. No.	5092



Tubes

	tubes ²⁾								tubes with screw cap								5127 ²⁴⁾
capacity in ml	5	7	15	25	50	100	100	250	15	50	12	25	30	25	25	250	
Ø x L in mm	13 x 75	12 x 100	17 x 100	24 x 100	34 x 100	44 x 100	40 x 115	65 x 115	17 x 120	29 x 115	17 x 100	25 x 90	25 x 110	25 x 90	61 x 122		
max. RCF ²⁾	3,005	3,005	2,952	2,898	2,952	2,952	2,952	3,095	3,095	3,095	3,005	2,826	2,898	2,826	3,095		
radius in mm	168	168	165	162	165	165	165	173	173	173	168	158	162	158	173		



Adapter

boring Ø x L in mm	12.8 x 42	13.2 x 82	17.5 x 56.7	25.5 x 74	35.5 x 77.5	45.5 x 76.5	42 x 76.5	66 x 103	17 x 85	30 x 85	17.5 x 79.5	26 x 73	25.5 x 74	26 x 73	62 x 90		
tubes per rotor	48	48	32	16	4	4	4	8	28	8	28	12	16	12	4		
Cat. No.	5128	5120	5136	5122	5124	5125	5126	1791	5129	5123	5121	5134	5122	5134	6319		

24)



Tubes

capacity in ml	290
Ø x L in mm	62 x 137
max. RCF ²⁾	3,095
radius in mm	173



Adapter

boring Ø x L in mm	62 x 90
tubes per rotor	4
Cat. No.	6319



2) Please note that the RCF values indicated refer only to rotor performance. The max. permissible RCF of tubes used should be verified with the individual manufacturers. The max. RCF for glass tubes annotated with footnote 2) is 4,000.

2.1) Please note that the RCF values indicated refer only to rotor performance. The max. permissible RCF of vessels used should be verified with the individual manufacturers. The max. RCF for ASTM Tubes annotated with footnote 2) is 700.

3.11) When using these tubes, bucket 5051 cannot be closed with lid 5053.

4) Please remove the spacer.

5) Tested by the TÜV in conformity with DIN EN 61010, section 2-020.

15) Also available with decanting aid. (Cat. No. 5247-91 or 5248-91).

24) At temperatures of over +40 °C and / or when not filled to capacity, bottles may warp during centrifugation.

Certifications / Registrations

Product certification:



Product registration:



Company certifications:



Dimensions

