



SUMMER DEALS HAIER CO₂ INCUBATOR & BIOSAFETY CABINETS

**SUMMER
PRICES**

Model	Humidity Reservoir	Capacity	Shelves (standard/maximum)	SUMMER DEAL PRICE
HCP-80	Max. 1.75 L/Min. 0.5 L	80 L	3 / 7	Check with your Hettich Benelux Account Manager
HCP-168	Max. 3.5 L/Min. 0.5 L	170 L	3 / 11	
HCP-258	Max. 5.5 L/Min. 0.5 L	258 L	3 / 13	

- Efficient microbial filter
- Active airflow control technology, to ensure stability & homogeneity
- Door heating function, eliminates microbial contamination
- 180°C sterilization on all surfaces
- Fast delivery (depending on the model)



Model	Remark	LCD Display or LED Touch screen	Int. dimensions W x D x H in mm	SUMMER DEAL PRICE
HR700-IIA2	Excl. base frame Incl. electric sash window	LCD Display	600 x 550 x 540	Check with your Hettich Benelux Account Manager
HR700-DJ	Base frame	n.a.	Height 700	
HR900-IIA2-S	Excl. base frame	LCD Display	936 x 620 x 635	
HR900-IIA2-S	Incl. base frame	LCD Display	936 x 620 x max.1345	
HR900-IIA2	Incl. base frame with gas and vacuum tap	LCD Display	920 x 620 x 650	
HR1200-IIA2	Incl. base frame with gas and vacuum tap	LCD Display	1220 x 620 x 650	
HR1200-IIA2-S	Incl. base frame	LCD Display	1310 x 620 x 635	
HR1500-IIA2	Incl. base frame with gas and vacuum tap	LCD Display	1520 x 620 x 650	
HR1200-IIA2-X	Incl. base frame	Touch screen	1230 x 600 x 655	
HR1500-IIA2-X	Incl. base frame	Touch screen	1530 x 660 x 655	
HR1800-IIA2-X	Incl. base frame	Touch screen	1830 x 600 x 655	

All models:

- 70% downflow/30% exhaust;
- class 2 safety;
- standard with an UV-light.

Protection for the operator, product and environment.

TERMS FOR THIS SUMMER DEAL

- The promotion runs from 21 July 2025 - 30 December 2025 (for orders placed during this period).
- During the promotional period, the promotion applies to the models mentioned in this flyer.
- The discount is not transferable to other persons or companies.
- The discount cannot be combined with any other discounts or price agreements.
- The discount does not apply to other accessories that may be purchased.
- For exact descriptions and applications, please contact your account manager.
- The discount cannot be exchanged for cash.
- The delivery terms and conditions of Hettich Benelux apply to this promotion.
- This promotion is only valid in the Netherlands at Hettich Benelux.
- Hettich Benelux reserves the right to change these promotional terms and conditions during the promotional period.

