



MoPlant (Double) - Dairy analysis

- Reliable and long life.
- Safe to operate due to exhaust and overpressure systems.
- All-in-one mobile unit.
- High degree of accuracy (better than 0.02%).

www.hettichbenelux.com

SPECIFICATIONS Double Moplant:	
Exterior dimensions D x W x H (mm)	1053 x 2383 x 1427 mm (with 2 side tables)
Exterior dimensions D x W x H (mm)	1053 x 1680 x 1427 mm (without side tables)
External / internal construction	Stainless steel ANSI 304
Components:	
Centrifuge	1 piece
Hot plate	2 pieces
Vacuum chambers	2 pieces
Cooling chambers	2 pieces
Centrifuge speed	710 rpm (50 Hz or 850 rpm 60 Hz) fixed
Hot plate temperature range	30 - 170 °C (0,1 °C adjustability)
Hot plate capacity	8 dishes per hot plate
Vacuum oven temperature range	50 - 170 °C (0,1 °C adjustability)
Ambient cooling chamber	15 - 30 °C
Digital controllers:	
Centrifuge	1 timer (adjustable)
Hot plate	2 temperatures (adjustable)
Vacuum chambers	2 temperatures (adjustable) / 2 timers (adjustable) / 2 vacuum read out
Cooling chambers	1 temperature / 2 timers (adjustable)
Over pressure technical chamber	1 pressure
TECHNICAL REQUIREMENTS:	
Electrical supply	400V 50/60 Hz, 25A
Main power supply	CEE socket 16A 400V*
Water supply	3/4"
Supply air	Ø 114
Air exhaust	Ø 114
TRANSPORT:	
Unit weight	± 490 kg
Unit weight incl. side tables	± 510 kg
Unit weight incl. packaging (wooden export crate)	± 665 kg
ARTICLE NUMBERS:	
MoPlant Double with side tables	8651.0010
MoPlant Double without side tables	8651.0015
Burettes station	9000.0224
Shaker for Mojonnier flasks	upon request

* Other voltage supplies are available on request.

