



MoPlant (Single) - Dairy analysis

- Reliable and long life.
- Safe to operate due to vapour exhaust and overpressure systems.
- All-in-one mobile unit.
- High degree of accuracy (better than 0.02%).

www.hettichbenelux.com

SPECIFICATIONS Single Moplant:	
Exterior dimensions D x W x H (mm)	1053 x 1939 x 1427 mm (with 2 side tables)
Exterior dimensions D x W x H (mm)	1053 x 1236 x 1427 mm (without side tables)
External / internal construction	Stainless steel ANSI 304
Components:	
Centrifuge	1 piece
Hot plate	1 piece
Vacuum chambers	1 piece
Cooling chambers	1 piece
Digital controllers:	
Centrifuge	1 timer (adjustable)
Hot plate	1 temperature (adjustable)
Vacuum chambers	1 temperature (adjustable) / 1 timer (adjustable) / 1 vacuum read out
Cooling chambers	1 temperature / 1 timer (adjustable)
Over pressure technical chamber	1 pressure
TECHNICAL REQUIREMENTS:	
Electrical supply	400V 50/60 Hz, 16A
Main power supply	CEE socket 16A 400V*
Water supply	3/4"
Supply air	Ø 114
Air exhaust	Ø 114
TRANSPORT:	
Unit weight	± 380 kg
Unit weight incl. side tables	± 400 kg
Unit weight incl. packaging (wooden export crate)	± 534 kg
ARTICLE NUMBERS:	
MoPlant Single with side tables	8650.0010
MoPlant Single without side tables	8650.0015
Burettes station	9000.0224
Shaker for Mojonnier flasks	upon request

* Other voltage supplies are available on request.

